

Description

PEGRx™ 2 is a crystallization reagent kit designed to evaluate an array of polymers of varying molecular weight in a medium ionic strength environment in the presence of additives, salts, volatile organics, and polyols versus a wide range of pH. Polymer reagents include Polyethylene glycols and Polyethylene glycol monomethyl ethers. The polymer molecular weight range between 200 and 20,000 is evaluated in a medium ionic strength formulation. Ten different buffers are used to span the range of pH between 3.5 and 9. The primary screen variables are polymer type, polymer molecular weight, pH and secondary reagents which include additives, salts, volatile organics, and polyols. PEGRx 2 is a straightforward, effective, and efficient screen for determining preliminary crystallization conditions in a medium ionic strength, polymeric reagent formulation in the pH range 3.5 to 9. The choice of reagents, reagent concentration, and pH were determined from public and proprietary databases. The formulation is biased with a focus on appropriate reagent concentrations as well as pH and polymer utilization, but is also balanced to ensure a fair sampling of novel reagents. PEGRx 2 formulations are unique and do not overlap with formulations found in other Hampton Research screens.

Sample Preparation

The macromolecular sample should be homogenous, as pure as is practically possible (>95%) and free of amorphous and particulate material. Remove amorphous material by centrifugation or micro-filtration prior to use (1, 2, 3).

The recommended sample concentration is 5 to 25 mg/ml in water. Initially, the sample should be free of any unnecessary additives in order to observe the effect of the PEGRx 2 variables. Ideally, the initial screen should be performed with a sample which has been dialyzed against water although ligands, ions, reducing agents, or other additives should be present as required by the sample for solubility, stability, or activity.


Performing The Screen

The following procedure describes the use of the PEGRx 2 with the Hanging Drop Vapor Diffusion method. The PEGRx 2 is also very compatible with the Sitting Drop, Sandwich Drop, MicroBatch, Microdialysis methods, and Free Interface Diffusion. A complete description of the Hanging, Sitting, Sandwich Drop, Dialysis and other crystallization methods are available from the Hampton Research Crystal Growth 101 Library.

1. Prepare a VDX Plate (HR3-140) for Hanging Drop Vapor Diffusion by applying a thin bead of cover slide sealant to the upper edge of each of the 24 reservoirs. One may also use a Greased VDX Plate (HR3-170). Forty-eight reservoirs are to be prepared for a complete PEGRx 2. See Figure 1.

Figure 1


Cross section of a reservoir in the VDX plate.



2. Using a clean pipet tip, pipet 1 ml of PEGRx 2 reagent 1 into reservoir A1. Discard the pipet tip, add a new pipet tip and pipet 1 ml of PEGRx 2 reagent 2 into reservoir A2. Repeat the procedure for the remaining 46 PEGRx 2 reagents using a clean pipet tip for each reagent so as to avoid reagent contamination and carry over.

3. Pipet 1 µl of the sample to the center of a clean, siliconized 22 mm diameter circle or square cover slide. See Figure 2.

Figure 2




4. Pipet 1 µl of PEGRx 2 reagent 1 from reservoir A1 into the sample droplet and mix by aspirating and dispensing the droplet several times, keeping the tip in the drop during mixing to avoid foaming. See Figure 2.

5. Working quickly to minimize evaporation, invert the cover slide and droplet over reservoir A1 and seal the cover slide onto the edge of the reservoir. See Figure 3.

Figure 3

Inverted siliconized coverslip placed over the reservoir.



6. Repeat operations 3 through 5 for the remaining 47 PEGRx 2 reagents.
7. If the quantity of sample permits, perform the PEGRx 2 in duplicate and incubate one set of plates at 4°C and the second set at room temperature. Incubate and store the crystallization plates in a stable temperature environment free of vibration.

Examine The Drop


Carefully examine the drops under a stereo microscope (10 to 100x magnification) immediately after setting up the screen. Record all observations and be particularly careful to scan the focal plane for small crystals. Observe the drops once each day for the first week, then once a week thereafter. Records should indicate whether the drop is clear, contains precipitate, and or crystals. It is helpful to describe the drop contents using descriptive terms. Adding magnitude is also helpful. Example: 4+ yellow/brown fine precipitate, 2+ small bipyramid crystals, clear drop, 3+ needle shaped crystals in 1+ white precipitate. One may also employ a standard numerical scoring scheme (Clear = 0, Precipitate = 1, Crystal = 10, etc). Figure 4 (on page 2) shows typical examples of what one might observe in a crystallization experiment.

Figure 4


Typical observations in a crystallization experiment




Clear Drop




Skin /
Precipitate




Precipitate




Precipitate /
Phase



Quasi
Crystals



Microcrystals




Needle
Cluster



Plates



Rod Cluster



Single
Crystal

Interpreting PEGRx 2

Clear drops indicate that either the relative supersaturation of the sample and reagent is too low or the drop has not yet completed equilibration. If the drop remains clear after 3 to 4 weeks consider repeating the PEGRx condition and doubling the sample concentration. If more than 33 of the 48 PEGRx 2 drops are clear consider doubling the sample concentration and repeating the entire screen.

Drops containing precipitate indicate that either the relative supersaturation of the sample and reagent is too high, the sample has denatured, or the sample is heterogeneous. To reduce the relative supersaturation, dilute the sample twofold and repeat the PEGRx 2 condition. If more than 33 of the 48 PEGRx 2 drops contain precipitate and no crystals are present, consider diluting the sample concentration in half and repeating the entire screen. If sample denaturation is suspect, take measures to stabilize the sample (add reducing agent, ligands, glycerol, salt, or other stabilizing agents). If the sample is impure, aggregated, or heterogeneous take measures to pursue homogeneity. It is possible to obtain crystals from precipitate so do not discard nor ignore a drop containing precipitate. If possible, examine drops containing precipitate under polarizing optics to differentiate precipitate from microcrystalline material.

If the drop contains a macromolecular crystal the relative supersaturation of the sample and reagent is good. The next step is to optimize the preliminary conditions (pH, polymer type, polymer concentration, polymer molecular weight, sample concentration, temperature, additives, and other crystallization variables) which produced the crystal in order to improve crystal size and quality.

Compare the observations between the 4°C and room temperature incubation to determine the effect of temperature on sample solubility. Different results in the same drops at different temperatures indicate that sample solubility is temperature dependent and that one should include temperature as a variable in subsequent screens and optimization experiments.

Retain and observe plates until the drops are dried out. Crystal growth can occur within 15 minutes or one year.

PEGRx 2 Formulation

PEGRx 2 reagents are formulated using the highest purity chemicals, ultrapure water (18.2 Megohm-cm, 5 ppb TOC) and are sterile filtered using 0.22 micron filters into sterile containers (no preservatives added).

PEGRx 2 reagents are readily reproduced using Hampton Research Optimize™ stock solutions of polymers and buffers. Optimize stock reagents make reproducing PEGRx 2 Screen reagents fast, convenient and easy. Dilutions can be performed directly into the crystallization plate using Optimize stock reagents.

No pH adjustments are made to PEGRx 2. Reagents are combined without further titration.

For further details about formulation, reproducing and optimizing reagents from the PEGRx 2 please refer to PEGRx 2 Fundamentals.

PEGRx 2 reagents are stable at room temperature and are best used before the "Best If Used By" date on the kit tubes. To enhance reagent stability it is recommended that PEGRx 2 be stored at 4°C or -20°C. Avoid ultraviolet light to preserve reagent stability.

References and Readings

1. Crystallization of Nucleic Acids and Proteins, Edited by A. Ducruix and R. Giegé, The Practical Approach Series, Oxford Univ. Press, 1992.
2. Current approaches to macromolecular crystallization. McPherson, A. Eur. J. Biochem. 189, 1-23, 1990.
3. Protein and Nucleic Acid Crystallization. Methods, A Companion to Methods in Enzymology, Academic Press, Volume 1, Number 1, August 1990.

Technical Support

Inquiries regarding PEGRx 2 reagent formulation, interpretation of screen results, optimization strategies and general inquiries regarding crystallization are welcome. Please e-mail, fax, or telephone your request to Hampton Research. Fax and e-mail Technical Support are available 24 hours a day. Telephone technical support is available 8:00 a.m. to 4:30 p.m. USA Pacific Standard Time.

Hampton Research
34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

Technical Support e-mail: tech@hrmail.com

Website: www.hamptonresearch.com

How to Reproduce PEGRx 2 Reagents

PEGRx™ 2 reagents and optimization conditions based on PEGRx 2 hits can be formulated using volumetric methods and carefully prepared reagent stocks (Table 1). Note the examples below.

Example 1. To prepare 1.0 milliliter of PEGRx 2 reagent 1 in a crystallization plate.

Solution Composition: 0.8 M Lithium sulfate monohydrate
0.1 M Sodium acetate trihydrate pH 4.0
4% v/v Polyethylene glycol 200

- 460 µl water³
- 400 µl 2.0 M Lithium sulfate monohydrate (CAS # 10102-25-7, Catalog # HR2-545)
- 100 µl 1.0 M Sodium acetate trihydrate pH 4.0 (CAS # 6131-90-4, Catalog # HR2-933-05)
- 40 µl 100% Polyethylene glycol 200 (CAS # 25322-68-3, Catalog # HR2-601)

Make no pH adjustments. Mix well by aspirating and dispensing the solution multiple times.

Example 2. To prepare 1.0 milliliter of PEGRx 2 reagent 15.

Solution Composition: 0.1 M Sodium malonate pH 8.0
0.1 M Tris pH 8.0
30% w/v Polyethylene glycol 1,000

- 570.6 µl water³
- 29.4 µl 3.4 M Sodium malonate pH 8.0 (CAS # 141-82-2, Catalog # HR2-807)
- 100 µl 1.0 M Tris pH 8.0 (CAS # 77-86-1, Catalog # HR2-900-11)
- 300 µl 100% Polyethylene glycol 1,000 (CAS # 25322-68-3, Catalog # HR2-523)

Make no pH adjustments. Mix well.

Example 3. To prepare 1 milliliter of PEGRx 2 reagent 48.

Solution Composition: 3% w/v Dextran sulfate sodium salt
0.1 M BICINE pH 8.5
15% w/v Polyethylene glycol 20,000

- 300 µl water³
- 100 µl 1.0 M BICINE pH 8.5 (CAS # 150-25-4, Catalog # HR2-999-10)

- 100 µl 30% w/v Dextran sulfate sodium salt (CAS # 9011-18-1, Catalog # HR2-428-51)
- 500 µl 30% w/v Polyethylene glycol 20,000 (CAS # 25322-68-3, Catalog # HR2-609)

Make no pH adjustments. Mix well.

³ ASTM Type I water.

Formulation Notes for PEGRx 2 Reagents

1. No additional pH adjustment is made to any reagent after formulation. Use the salts and buffers in Table 1 to reproduce a PEGRx 2 reagent.
2. All Optimize solutions and screen reagents are sterile filtered using 0.22 µm filters into sterile containers.
3. Add water first as this will help maintain the solubility of subsequently added reagents.
4. When formulating reagents using a pipet, add the largest volume last (except water). Use this larger volume setting to aspirate and dispense the reagent until the solution is mixed.
5. When formulating reagents using a pipet, use a clean, sterile pipet tip for each reagent added to the solution.
6. Use Optimize™, Custom Shop™, StockOptions™ pH, and StockOptions™ pH buffer kits from Hampton Research to systematically vary the pH as a crystallization variable.
7. The measured final pH of all PEGRx 2 reagents is available at www.hamptonresearch.com. Search using catalog number HR2-084 and follow the link to the 'PEGRx 2 pH and Conductivity' document.

pH as a Crystallization Variable

The buffers listed in Table 2 can be used to vary the pH as a crystallization variable and are recommended when optimizing a crystal grown from a PEGRx 2 kit.

Optimize™ buffer stocks are supplied as a 100 milliliters sterile filtered solution.

StockOptions™ buffer kits contain 10 milliliters each of ready to pipet buffers, titrated in 0.1 pH increments over the indicated pH range. The number of reagents offered in a StockOptions buffer kit depends upon the pH range of the buffer. The broader the pH range, the more buffers in the kit.

Custom Shop™ ready to pipet buffers are made to order, pH titrated buffer stocks from Hampton Research.

Online Information

Visit www.hamptonresearch.com and enter one of the following:

- Reagent Catalog Number
- Kit Catalog Number
- CAS Number
- Reagent Name

To obtain reagent specifications, pH titration tables, user guides, certificates of analysis, material safety data sheets (MSDS), and any other additional information.

MakeTray™

MakeTray is a free, web based program at www.hamptonresearch.com which generates both a pipetting worksheet and a reagent formulation document for crystallization set ups. MakeTray allows one to enter general information about the sample and experiment, which is then printed on the pipet worksheet and the reagent formulation document. The plate size can be customized for any number of wells, so MakeTray works for: 24, 48, and 96 well plates. MakeTray is especially useful for the design and formulation of crystal optimization experiments.

Table 1. Recommended reagents for the formulation of PEGRx 2 and optimization reagents.

Each of these reagents are available as an Optimize™ crystallization grade reagent from Hampton Research. Table 1 provides the common chemical name, the Hampton Research catalog number, supplied stock concentration, the supplied volume, and the CAS number for each reagent. For more information on a specific Optimize reagent, go to

www.hamptonresearch.com. Using Search, enter either the catalog number, CAS number, or chemical name to obtain additional information for the Optimize reagent, including a Certificate of Analysis and MSDS (where applicable).

| Additive | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
|--------------------------------|----------------------------|--------------------|-----------------|-------------|
| CYMAL®-7 | HR2-406-23 | 0.15 mM | 1.0 ml | 349477-49-2 |
| Dextran sulfate sodium salt | HR2-428-51 | 30% w/v | 1.0 ml | 9011-18-1 |
| n-Octyl-β-D-glucoside | HR2-428-71 | 5% w/v | 1.0 ml | 29836-26-8 |
| | | | | |
| Organics (Volatile) | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
| 1,4-Dioxane | HR2-617 | 100% | 200 ml | 123-91-1 |
| 2-Propanol | HR2-619 | 100% | 200 ml | 67-63-0 |
| | | | | |
| Polyol | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
| (+/-)-2-Methyl-2,4-pentanediol | HR2-627 | 100% | 200 ml | 107-41-5 |
| Ethylene glycol | HR2-621 | 100% | 100 ml | 107-21-1 |
| Polyethylene glycol 200 | HR2-601 | 100% | 200 ml | 25322-68-3 |
| Polyethylene glycol 400 | HR2-603 | 100% | 200 ml | 25322-68-3 |
| | | | | |
| Salt | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
| Ammonium acetate | HR2-565 | 1.0 M | 100 ml | 631-61-8 |
| | HR2-799 | 8.0 M | 200 ml | 631-61-8 |
| | | | | |
| (Salt Continued on page 3) | | | | |

Table 1 (Continued). Recommended reagents for the formulation of PEGRx 2 and optimization reagents.

| Salt | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
|--|----------------------------|--------------------|-----------------|-------------|
| Ammonium citrate tribasic pH 7.0 | HR2-759 | 2.5 M | 200 ml | 3458-72-8 |
| Ammonium sulfate | HR2-541 | 3.5 M | 200 ml | 7783-20-2 |
| Cadmium chloride hydrate | HR2-715 | 1.0 M | 100 ml | 654054-66-7 |
| Calcium chloride dihydrate | HR2-557 | 2.0 M | 100 ml | 10035-04-8 |
| L-Proline | HR2-775 | 1.0 M | 100 ml | 147-85-3 |
| Lithium sulfate monohydrate | HR2-545 | 2.0 M | 200 ml | 10377-48-7 |
| Magnesium chloride hexahydrate | HR2-559 | 2.0 M | 100 ml | 7791-18-6 |
| | HR2-803 | 5.0 M | 200 ml | 7791-18-6 |
| Magnesium formate dihydrate | HR2-537 | 1.0 M | 200 ml | 557-39-1 |
| DL-Malic acid pH 7.0 | HR2-761 | 3.0 M | 200 ml | 6915-15-7 |
| Nickel(II) chloride hexahydrate | HR2-687 | 4.0 M | 200 ml | 7791-20-0 |
| Potassium formate | HR2-667 | 14.0 M | 200 ml | 590-29-4 |
| Potassium sodium tartrate tetrahydrate | HR2-539 | 1.5 M | 200 ml | 6381-59-5 |
| Succinic acid pH 7.0 | HR2-709 | 1.2 M | 200 ml | 110-15-6 |
| Sodium chloride | HR2-637 | 5.0 M | 200 ml | 7647-14-5 |
| Sodium formate | HR2-547 | 7.0 M | 200 ml | 141-53-7 |
| Sodium malonate pH 5.0 | HR2-749 | 3.4 M | 200 ml | 141-82-2 |
| Sodium malonate pH 6.0 | HR2-751 | 3.4 M | 200 ml | 141-82-2 |
| Sodium malonate pH 8.0 | HR2-807 | 3.4 M | 200 ml | 141-82-2 |
| Taccimate pH 4.0 | HR2-823 | 100% | 200 ml | N/A |
| Taccimate pH 6.0 | HR2-827 | 100% | 200 ml | N/A |
| Taccimate pH 7.0 | HR2-755 | 100% | 200 ml | N/A |
| | | | | |
| Polymer | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
| Polyethylene glycol 200 | HR2-601 | 100% | 200 ml | 25322-68-3 |
| Polyethylene glycol 400 | HR2-603 | 100% | 200 ml | 25322-68-3 |
| Polyethylene glycol 1,000 | HR2-523 | 50% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol 1,500 | HR2-525 | 50% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol 3,350 | HR2-527 | 50% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol 4,000 | HR2-529 | 50% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol 6,000 | HR2-533 | 50% w/v | 200 ml | 25322-68-3 |
| | | | | |

(Polymer Continued on page 4)

Table 1 (Continued). Recommended reagents for the formulation of PEGRx 2 and optimization reagents.

| Polymer | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
|--|----------------------------|--------------------|-----------------|-------------|
| Polyethylene glycol 8,000 | HR2-535 | 50% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol 10,000 | HR2-607 | 50% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol 20,000 | HR2-609 | 30% w/v | 200 ml | 25322-68-3 |
| Polyethylene glycol monomethyl ether 550 | HR2-611 | 100% | 200 ml | 9004-74-4 |
| Polyethylene glycol monomethyl ether 2,000 | HR2-613 | 50% w/v | 200 ml | 9004-74-4 |
| Polyethylene glycol monomethyl ether 5,000 | HR2-615 | 50% w/v | 200 ml | 9004-74-4 |
| | | | | |
| Buffer | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # |
| BICINE pH 8.5 ² | HR2-999-10 | 1.0 M | 185 ml | 150-25-4 |
| BIS-TRIS pH 6.5 ¹ | HR2-783 | 1.0 M | 100 ml | 6976-37-0 |
| BIS-TRIS propane pH 9.0 ¹ | HR2-993-28 | 1.0 M | 185 ml | 64431-96-5 |
| Citric acid pH 3.5 ² | HR2-757 | 1.0 M | 100 ml | 77-92-9 |
| HEPES pH 7.5 ² | HR2-729 | 1.0 M | 100 ml | 7365-45-9 |
| Imidazole pH 7.0 ¹ | HR2-819 | 1.0 M | 100 ml | 288-32-4 |
| MES monohydrate pH 6.0 ² | HR2-943-09 | 1.0 M | 185 ml | 145224-94-8 |
| Sodium acetate trihydrate pH 4.0 ¹ | HR2-933-05 | 1.0 M | 185 ml | 6131-90-4 |
| Sodium acetate trihydrate pH 4.5 ¹ | HR2-789 | 1.0 M | 100 ml | 6131-90-4 |
| Sodium citrate tribasic dihydrate pH 5.0 ¹ | HR2-935-09 | 1.0 M | 185 ml | 6132-04-3 |
| Sodium citrate tribasic dihydrate pH 5.5 ¹ | HR2-935-14 | 1.0 M | 185 ml | 6132-04-3 |
| Tris pH 8.0 ¹ | HR2-900-11 | 1.0 M | 185 ml | 77-86-1 |
| | | | | |
| ¹ pH titrated using Hydrochloric acid (HR2-581) CAS # 7647-01-0 | | | | |
| ² pH titrated using Sodium hydroxide (HR2-583) CAS # 1310-73-2 | | | | |

Table 2. Recommended buffers for screening the pH of PEGRx 2 and optimization reagents.

| Buffer Solution or Kit | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # | pH range |
|---|----------------------------|--------------------|-----------------|------------|-----------|
| StockOptions™ BICINE | HR2-999-*** | 1.0 M | 185 ml | 150-25-4 | 7.6 - 9.0 |
| StockOptions™ Bis-Tris kit ⁴ | HR2-106 | 1.0 M | 10 ml each | 6976-37-0 | 5.5 - 7.5 |
| StockOptions™ Bis-Tris propane | HR2-993-*** | 1.0 M | 185 ml | 64431-96-5 | 6.3 - 9.5 |
| | | | | | |

(Buffer Solution or Kit Continued on page 5)

Table 2. Recommended buffers for screening the pH of PEGRx 2 and optimization reagents.

| Buffer Solution or Kit | Hampton Research Catalog # | Supplied [Stock] | Supplied Volume | CAS # | pH range |
|---|----------------------------|--------------------|-----------------|-------------|-----------|
| StockOptions™ Citric acid ⁴ | HR2-104 | 1.0 M | 10 ml each | 77-92-9 | 2.2 - 6.5 |
| StockOptions™ Hepes kit ⁴ | HR2-102 | 1.0 M | 10 ml each | 7365-45-9 | 6.8 - 8.2 |
| StockOptions™ Imidazole | HR2-995-** | 1.0 M | 185 ml | 288-32-4 | 6.2 - 7.8 |
| StockOptions™ MES kit ⁴ | HR2-243 | 1.0 M | 10 ml each | 145224-94-8 | 5.2 - 7.1 |
| StockOptions™ Sodium Acetate kit ⁴ | HR2-233 | 1.0 M | 10 ml each | 6131-90-4 | 3.6 - 5.6 |
| StockOptions™ Sodium Citrate kit ⁴ | HR2-235 | 1.0 M | 10 ml each | 6132-04-3 | 4.2 - 6.5 |
| StockOptions™ Tris ⁴ | HR2-100 | 1.0 M | 10 ml each | 77-86-1 | 7.0 - 9.0 |
| | | | | | |

⁴ Individual StockOptions buffers titrated to any pH within the kit's pH range are available in 185 ml volumes from the Hampton Research Custom Shop™

** Refers to the reagent number in the kit. For example, reagent number 1 = HR2-993-01 (pH 6.3)

Technical Support

Inquiries regarding PEGRx 2 Fundamentals, interpretation of screen results, optimization strategies and general inquiries regarding crystallization are welcome. Please e-mail, fax, or telephone your request to Hampton Research. Fax and e-mail Technical Support are available 24 hours a day. Telephone technical support is available 8:00 a.m. to 4:30 p.m. USA Pacific Standard Time.

Hampton Research
34 Journey
Aliso Viejo, CA 92656-3317 U.S.A.
Tel: (949) 425-1321 • Fax: (949) 425-1611
Technical Support e-mail: tech@hrmail.com
Website: www.hamptonresearch.com

| Tube # | Additive / Salt / Volatile Organic / Polyol | Tube # | Buffer ◇ | Tube # | Polymer |
|--------|--|--------|--|--------|--|
| 1. | 0.8 M Lithium sulfate monohydrate | 1. | 0.1 M Sodium acetate trihydrate pH 4.0 | 1. | 4% v/v Polyethylene glycol 200 |
| 2. | 0.2 M Lithium sulfate monohydrate | 2. | 0.1 M Sodium citrate tribasic dihydrate pH 5.0 | 2. | 26% v/v Polyethylene glycol 200 |
| 3. | 0.05 M Calcium chloride dihydrate | 3. | 0.1 M MES monohydrate pH 6.0 | 3. | 45% v/v Polyethylene glycol 200 |
| 4. | 28% v/v 2-Propanol | 4. | 0.1 M BIS-TRIS pH 6.5 | 4. | 3% v/v Polyethylene glycol 200 |
| 5. | 20% v/v Tacsimate pH 7.0 | 5. | 0.1 M HEPES pH 7.5 | 5. | 2% v/v Polyethylene glycol 200 |
| 6. | 10% v/v 2-Propanol | 6. | 0.1 M Sodium citrate tribasic dihydrate pH 5.0 | 6. | 26% v/v Polyethylene glycol 400 |
| 7. | 0.2 M Ammonium acetate | 7. | 0.1 M Sodium citrate tribasic dihydrate pH 5.5 | 7. | 24% v/v Polyethylene glycol 400 |
| 8. | 0.2 M Ammonium sulfate | 8. | 0.1 M BIS-TRIS pH 6.5 | 8. | 18% v/v Polyethylene glycol 400 |
| 9. | 0.19 mM CYMAL®-7 | 9. | 0.1 M HEPES pH 7.5 | 9. | 40% v/v Polyethylene glycol 400 |
| 10. | 6% v/v 2-Propanol | 10. | 0.1 M Sodium acetate trihydrate pH 4.5 | 10. | 26% v/v Polyethylene glycol monomethyl ether 550 |
| 11. | 1.8 M Ammonium sulfate | 11. | 0.1 M BIS-TRIS pH 6.5 | 11. | 2% v/v Polyethylene glycol monomethyl ether 550 |
| 12. | 0.15 M DL-Malic acid pH 7.0 | 12. | 0.1 M Imidazole pH 7.0 | 12. | 22% v/v Polyethylene glycol monomethyl ether 550 |
| 13. | 0.1 M Succinic acid pH 7.0 | 13. | 0.1 M BICINE pH 8.5 | 13. | 30% v/v Polyethylene glycol monomethyl ether 550 |
| 14. | 0.1 M Lithium sulfate monohydrate | 14. | 0.1 M Sodium citrate tribasic dihydrate pH 5.5 | 14. | 20% w/v Polyethylene glycol 1,000 |
| 15. | 0.1 M Sodium malonate pH 8.0 | 15. | 0.1 M Tris pH 8.0 | 15. | 30% w/v Polyethylene glycol 1,000 |
| 16. | 4% v/v (+/-)-2-Methyl-2,4-pentanediol | 16. | 0.1 M Citric acid pH 3.5 | 16. | 20% w/v Polyethylene glycol 1,500 |
| 17. | 0.2 M L-Proline | 17. | 0.1 M HEPES pH 7.5 | 17. | 24% w/v Polyethylene glycol 1,500 |
| 18. | 10% v/v 2-Propanol | 18. | 0.1 M BICINE pH 8.5 | 18. | 30% w/v Polyethylene glycol 1,500 |
| 19. | 0.1 M Sodium chloride | 19. | 0.1 M BIS-TRIS propane pH 9.0 | 19. | 25% w/v Polyethylene glycol 1,500 |
| 20. | 0.02 M Nickel(II) chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate, 0.02 M Cadmium chloride hydrate | 20. | 0.1 M Sodium acetate trihydrate pH 4.5 | 20. | 24% w/v Polyethylene glycol monomethyl ether 2,000 |
| 21. | 20% v/v 2-Propanol | 21. | 0.1 M MES monohydrate pH 6.0 | 21. | 20% w/v Polyethylene glycol monomethyl ether 2,000 |
| 22. | 0.2 M Ammonium citrate tribasic pH 7.0 | 22. | 0.1 M Imidazole pH 7.0 | 22. | 20% w/v Polyethylene glycol monomethyl ether 2,000 |
| 23. | 4.0 M Potassium formate | 23. | 0.1 M BIS-TRIS propane pH 9.0 | 23. | 2% w/v Polyethylene glycol monomethyl ether 2,000 |
| 24. | 50% v/v Tacsimate pH 4.0 | 24. | 0.1 M Sodium acetate trihydrate pH 4.5 | 24. | 1% w/v Polyethylene glycol 3,350 |
| 25. | 0.10% w/v n-Octyl-β-D-glucoside | 25. | 0.1 M Sodium citrate tribasic dihydrate pH 5.5 | 25. | 22% w/v Polyethylene glycol 3,350 |
| 26. | 2% v/v Tacsimate pH 7.0, 5% v/v 2-Propanol | 26. | 0.1 M Imidazole pH 7.0 | 26. | 8% w/v Polyethylene glycol 3,350 |
| 27. | 2% v/v 1,4-Dioxane | 27. | 0.1 M Tris pH 8.0 | 27. | 15% w/v Polyethylene glycol 3,350 |
| 28. | 18% v/v 2-Propanol | 28. | 0.1 M Sodium citrate tribasic dihydrate pH 5.5 | 28. | 20% w/v Polyethylene glycol 4,000 |
| 29. | 6% v/v Tacsimate pH 6.0 | 29. | 0.1 M MES monohydrate pH 6.0 | 29. | 25% w/v Polyethylene glycol 4,000 |
| 30. | 0.2 M Magnesium formate dihydrate | 30. | 0.1 M Sodium acetate trihydrate pH 4.0 | 30. | 18% w/v Polyethylene glycol monomethyl ether 5,000 |
| 31. | 2% v/v Polyethylene glycol 400 | 31. | 0.1 M Imidazole pH 7.0 | 31. | 24% w/v Polyethylene glycol monomethyl ether 5,000 |
| 32. | 0.2 M Sodium formate | 32. | 0.1 M BICINE pH 8.5 | 32. | 20% w/v Polyethylene glycol monomethyl ether 5,000 |
| 33. | 4% v/v 2-Propanol | 33. | 0.1 M BIS-TRIS propane pH 9.0 | 33. | 20% w/v Polyethylene glycol monomethyl ether 5,000 |
| 34. | 6% v/v Ethylene glycol | 34. | 0.1 M Citric acid pH 3.5 | 34. | 10% w/v Polyethylene glycol 6,000 |
| 35. | 0.15 M Lithium sulfate monohydrate | 35. | 0.1 M Citric acid pH 3.5 | 35. | 18% w/v Polyethylene glycol 6,000 |
| 36. | 10% v/v 2-Propanol | 36. | 0.1 M Sodium acetate trihydrate pH 4.0 | 36. | 22% w/v Polyethylene glycol 6,000 |
| 37. | 0.2 M Sodium chloride | 37. | 0.1 M Sodium acetate trihydrate pH 4.0 | 37. | 22% w/v Polyethylene glycol 8,000 |
| 38. | 20% v/v 2-Propanol | 38. | 0.1 M Tris pH 8.0 | 38. | 5% w/v Polyethylene glycol 8,000 |
| 39. | 10% v/v Polyethylene glycol 200 | 39. | 0.1 M BIS-TRIS propane pH 9.0 | 39. | 18% w/v Polyethylene glycol 8,000 |
| 40. | 15% v/v 2-Propanol | 40. | 0.1 M Sodium citrate tribasic dihydrate pH 5.0 | 40. | 10% w/v Polyethylene glycol 10,000 |
| 41. | 0.4 M Sodium malonate pH 6.0 | 41. | 0.1 M MES monohydrate pH 6.0 | 41. | 0.5% w/v Polyethylene glycol 10,000 |
| 42. | 0.2 M Potassium sodium tartrate tetrahydrate | 42. | 0.1 M BIS-TRIS pH 6.5 | 42. | 10% w/v Polyethylene glycol 10,000 |
| 43. | 5% v/v (+/-)-2-Methyl-2,4-pentanediol | 43. | 0.1 M HEPES pH 7.5 | 43. | 10% w/v Polyethylene glycol 10,000 |
| 44. | 0.2 M Ammonium acetate | 44. | 0.1 M Tris pH 8.0 | 44. | 16% w/v Polyethylene glycol 10,000 |
| 45. | 5% v/v 2-Propanol | 45. | 0.1 M Citric acid pH 3.5 | 45. | 6% w/v Polyethylene glycol 20,000 |
| 46. | 1.0 M Sodium malonate pH 5.0 | 46. | 0.1 M Sodium acetate trihydrate pH 4.5 | 46. | 2% w/v Polyethylene glycol 20,000 |
| 47. | 0.2 M Magnesium chloride hexahydrate | 47. | 0.1 M Sodium citrate tribasic dihydrate pH 5.0 | 47. | 10% w/v Polyethylene glycol 20,000 |
| 48. | 3% w/v Dextran sulfate sodium salt | 48. | 0.1 M BICINE pH 8.5 | 48. | 15% w/v Polyethylene glycol 20,000 |

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

PEGRx™ 2 contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.

Sample: _____ Sample Concentration: _____
 Sample Buffer: _____ Date: _____
 Reservoir Volume: _____ Temperature: _____
 Drop Volume: Total _____ µl Sample _____ µl Reservoir _____ µl Additive _____ µl

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals

- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2 mm)
- 9 Single Crystals (3D Growth > 0.2 mm)

PEGRx™ 2 - HR2-084 Scoring Sheet

Date: Date: Date:

| | | | | |
|-----|--|--|--|--|
| 1. | 0.8 M Lithium sulfate monohydrate, 0.1 M Sodium acetate trihydrate pH 4.0, 4% v/v Polyethylene glycol 200 | | | |
| 2. | 0.2 M Lithium sulfate monohydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.0, 26% v/v Polyethylene glycol 200 | | | |
| 3. | 0.05 M Calcium chloride dihydrate, 0.1 M MES monohydrate pH 6.0, 45% v/v Polyethylene glycol 200 | | | |
| 4. | 28% v/v 2-Propanol, 0.1 M BIS-TRIS pH 6.5, 3% v/v Polyethylene glycol 200 | | | |
| 5. | 20% v/v Tacsimate pH 7.0, 0.1 M HEPES pH 7.5, 2% v/v Polyethylene glycol 200 | | | |
| 6. | 10% v/v 2-Propanol, 0.1 M Sodium citrate tribasic dihydrate pH 5.0, 26% v/v Polyethylene glycol 400 | | | |
| 7. | 0.2 M Ammonium acetate, 0.1 M Sodium citrate tribasic dihydrate pH 5.5, 24% v/v Polyethylene glycol 400 | | | |
| 8. | 0.2 M Ammonium sulfate, 0.1 M BIS-TRIS pH 6.5, 18% v/v Polyethylene glycol 400 | | | |
| 9. | 0.19 mM CYMAL®-7, 0.1 M HEPES pH 7.5, 40% v/v Polyethylene glycol 400 | | | |
| 10. | 6% v/v 2-Propanol, 0.1 M Sodium acetate trihydrate pH 4.5, 26% v/v Polyethylene glycol monomethyl ether 550 | | | |
| 11. | 1.8 M Ammonium sulfate, 0.1 M BIS-TRIS pH 6.5, 2% v/v Polyethylene glycol monomethyl ether 550 | | | |
| 12. | 0.15 M DL-Malic acid pH 7.0, 0.1 M Imidazole pH 7.0, 22% v/v Polyethylene glycol monomethyl ether 550 | | | |
| 13. | 0.1 M Succinic acid pH 7.0, 0.1 M BICINE pH 8.5, 30% v/v Polyethylene glycol monomethyl ether 550 | | | |
| 14. | 0.1 M Lithium sulfate monohydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.5, 20% w/v Polyethylene glycol 1,000 | | | |
| 15. | 0.1 M Sodium malonate pH 8.0, 0.1 M Tris pH 8.0, 30% w/v Polyethylene glycol 1,000 | | | |
| 16. | 4% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.1 M Citric acid pH 3.5, 20% w/v Polyethylene glycol 1,500 | | | |
| 17. | 0.2 M L-Proline, 0.1 M HEPES pH 7.5, 24% w/v Polyethylene glycol 1,500 | | | |
| 18. | 10% v/v 2-Propanol, 0.1 M BICINE pH 8.5, 30% w/v Polyethylene glycol 1,500 | | | |
| 19. | 0.1 M Sodium chloride, 0.1 M BIS-TRIS propane pH 9.0, 25% w/v Polyethylene glycol 1,500 | | | |
| 20. | 0.02 M Nickel(II) chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate, 0.02 M Cadmium chloride hydrate, 0.1 M Sodium acetate trihydrate pH 4.5, 24% w/v Polyethylene glycol monomethyl ether 2,000 | | | |
| 21. | 20% v/v 2-Propanol, 0.1 M MES monohydrate pH 6.0, 20% w/v Polyethylene glycol monomethyl ether 2,000 | | | |
| 22. | 0.2 M Ammonium citrate tribasic pH 7.0, 0.1 M Imidazole pH 7.0, 20% w/v Polyethylene glycol monomethyl ether 2,000 | | | |
| 23. | 4.0 M Potassium formate, 0.1 M BIS-TRIS propane pH 9.0, 2% w/v Polyethylene glycol monomethyl ether 2,000 | | | |
| 24. | 50% v/v Tacsimate pH 4.0, 0.1 M Sodium acetate trihydrate pH 4.5, 1% w/v Polyethylene glycol 3,350 | | | |
| 25. | 0.10% w/v n-Octyl-β-D-glucoside, 0.1 M Sodium citrate tribasic dihydrate pH 5.5, 22% w/v Polyethylene glycol 3,350 | | | |
| 26. | 2% v/v Tacsimate pH 7.0, 5% v/v 2-Propanol, 0.1 M Imidazole pH 7.0, 8% w/v Polyethylene glycol 3,350 | | | |
| 27. | 2% v/v 1,4-Dioxane, 0.1 M Tris pH 8.0, 15% w/v Polyethylene glycol 3,350 | | | |
| 28. | 18% v/v 2-Propanol, 0.1 M Sodium citrate tribasic dihydrate pH 5.5, 20% w/v Polyethylene glycol 4,000 | | | |
| 29. | 6% v/v Tacsimate pH 6.0, 0.1 M MES monohydrate pH 6.0, 25% w/v Polyethylene glycol 4,000 | | | |
| 30. | 0.2 M Magnesium formate dihydrate, 0.1 M Sodium acetate trihydrate pH 4.0, 18% w/v Polyethylene glycol monomethyl ether 5,000 | | | |
| 31. | 2% v/v Polyethylene glycol 400, 0.1 M Imidazole pH 7.0, 24% w/v Polyethylene glycol monomethyl ether 5,000 | | | |
| 32. | 0.2 M Sodium formate, 0.1 M BICINE pH 8.5, 20% w/v Polyethylene glycol monomethyl ether 5,000 | | | |
| 33. | 4% v/v 2-Propanol, 0.1 M BIS-TRIS propane pH 9.0, 20% w/v Polyethylene glycol monomethyl ether 5,000 | | | |
| 34. | 6% v/v Ethylene glycol, 0.1 M Citric acid pH 3.5, 10% w/v Polyethylene glycol 6,000 | | | |
| 35. | 0.15 M Lithium sulfate monohydrate, 0.1 M Citric acid pH 3.5, 18% w/v Polyethylene glycol 6,000 | | | |
| 36. | 10% v/v 2-Propanol, 0.1 M Sodium acetate trihydrate pH 4.0, 22% w/v Polyethylene glycol 6,000 | | | |
| 37. | 0.2 M Sodium chloride, 0.1 M Sodium acetate trihydrate pH 4.0, 22% w/v Polyethylene glycol 8,000 | | | |
| 38. | 20% v/v 2-Propanol, 0.1 M Tris pH 8.0, 5% w/v Polyethylene glycol 8,000 | | | |
| 39. | 10% v/v Polyethylene glycol 200, 0.1 M BIS-TRIS propane pH 9.0, 18% w/v Polyethylene glycol 8,000 | | | |
| 40. | 15% v/v 2-Propanol, 0.1 M Sodium citrate tribasic dihydrate pH 5.0, 10% w/v Polyethylene glycol 10,000 | | | |
| 41. | 0.4 M Sodium malonate pH 6.0, 0.1 M MES monohydrate pH 6.0, 0.5% w/v Polyethylene glycol 10,000 | | | |
| 42. | 0.2 M Potassium sodium tartrate tetrahydrate, 0.1 M BIS-TRIS pH 6.5, 10% w/v Polyethylene glycol 10,000 | | | |
| 43. | 5% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.1 M HEPES pH 7.5, 10% w/v Polyethylene glycol 10,000 | | | |
| 44. | 0.2 M Ammonium acetate, 0.1 M Tris pH 8.0, 16% w/v Polyethylene glycol 10,000 | | | |
| 45. | 5% v/v 2-Propanol, 0.1 M Citric acid pH 3.5, 6% w/v Polyethylene glycol 20,000 | | | |
| 46. | 1.0 M Sodium malonate pH 5.0, 0.1 M Sodium acetate trihydrate pH 4.5, 2% w/v Polyethylene glycol 20,000 | | | |
| 47. | 0.2 M Magnesium chloride hexahydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.0, 10% w/v Polyethylene glycol 20,000 | | | |
| 48. | 3% w/v Dextran sulfate sodium salt, 0.1 M BICINE pH 8.5, 15% w/v Polyethylene glycol 20,000 | | | |