

**Crystallization Screening using
Trace Fluorescent Labeling**


Or

**How Using Fluorescence can be a Brilliant
Answer to YOUR Crystallization Problem.**


Marc L. Pusey
iXpressGenes, Inc.

Even Heavily Labeled Proteins Crystallize

Original impetus was to produce FRET pairs for use in studying the lysozyme crystal nucleation process. Lysozyme was site-specifically labeled on asp-101, then Cation-exchange purified to remove labeled from unlabeled protein. Crystallization experiments were carried out to determine if the labeled proteins would indeed crystallize.



Cascade blue-lys
Ex = 396 nm, Em = 410 nm




Lucifer yellow – lys
Ex = 428 nm, Em = 536 nm



Texas red - UCAN

The proteins in the above images were > 50% (cascade blue and lucifer yellow-lys to > 90%) labeled with the indicated fluorescent probe. This gave an early indication that protein could be labeled and would still crystallize. The fluorescent microscope was set up for green fluorescent protein.

Evolution of the Trace Fluorescence Idea - Implementation



Frustration led to development of an in-house assembled scanning microscope, that imaged every well and saved them to disk. Supported by NIH Phase I STTR.



Original work used an Olympus Stereo Fluorescence microscope. Had to manually grab and save each image. Supported by NIH R21 grant while I was at NASA.

The in-house design has been developed by a local design engineering company to a more presentable product. Currently supports 2 fluorescent and white light imaging, or 3 fluorescent and no white light imaging. ~15 minutes to scan a 96 well plate, 3 drops/well, single wavelength.



Fluorescence Methods

Different Approaches

Native tryptophan fluorescence

Positives – protein is ‘untouched’

Negatives – requires UV optics, illumination, therefore much higher cost

- cannot use a ‘hands (eyes) on’ approach.**
- not all proteins have a tryptophan.**

Adding free probe to crystallizing solution

Positives – Easy to implement

- no covalent modification needed**

Negatives – possible uncertainty about where probe is partitioning

- cannot be used with IMPs**
- less sensitive**

Covalent modification of protein with probe

Positives – Certainty about following the protein

Negatives – Most work to use

- Worries about effects on nucleation, crystal quality, inhibit use.**

Other benefits found using the covalent modification approach

That may or may not be found with the other two.

Trace Fluorescent Labeling

The methodology to be described is TRACE fluorescent labeling, where on average **ONLY ~ 0.1 – 0.2 %** of the protein molecules have a fluorescent probe covalently attached. This is 1 or 2 molecules per thousand being modified. IF your protein was 100% pure (it almost certainly isn't) then this would drop the “purity” to 99.9 or 99.8 %.

The basis of the technique is:

Fluorescence intensity is proportional to fluor concentration.

The highest solute concentration is in the crystalline state.


Therefore crystals are brighter objects.

This gives the guiding rule → **INTENSITY = STRUCTURE**

Note the emphasis on **TRACE**, meaning < 0.5% (and preferably Lower. NOT 1%, 5%, 10%, 50%, or 100%, but < 0.5%!!)

Sensitivity – How Low?

Glucose Isomerase and Xylanase were derivatized with CR, then dilutions made in 1/5x steps and crystallization plates set up. For a 2nd level of comparison, Crystals were also set up using $2 \times 10^{-5} \text{M}$ ANS (8-Anilino-1-naphthalenesulfonic acid) using the IZIT approach.



The intensity of the excitation source also affects the lower detection limit. The above images were collected using a single LED for the source illumination. By now using a higher powered LED source I estimate the sensitivity can be adjusted ~ an order of magnitude lower. A more sensitive camera would lower this still more.

Trace Labeling is Simple and Quick

A procedure to quickly prepare trace fluorescently labeled protein solution.

Start with -
0.8 mL of protein
@ 125% desired
final conc.



Remove aliquot of
0.08 mL, buffer
exchange if
necessary
~2 min.



Typically takes 5 to 30 minutes and
is found to yield ~0.2 % derivatized
protein.

Add 0.5-2.0 μ L of fluorescent
probe solution (more with
higher protein MWs)

Add back to the
Starting solution,
Bring volume to 1.0 ml
DONE!!

Buffer exchange
(to stock solution
buffer) ~2 min.




Incubate while
reequilibrating
the desalting
column (if using
a single desalting
column) ~10 min.

This procedure can be
scaled up or down as
needed. To reduce the
amount of protein the
reactive probe solution
should be diluted. Smaller
centrifugal desalting columns
are available.

WHY DO IT?

If intensity is proportional to structure – find the crystals & conditions of interest below.



The BIG QUESTION

Does the presence of the fluorescent probe affect the crystallization?

The short answer is – **NO!**

(but that won't stop me from talking about it in detail for the next couple of minutes)

Experiments to again¹ determine the answer (now in progress)

Using a series of test proteins, some previously worked on, some not:

Using a single protein preparation –

Prepare labeled and unlabeled protein at same concentrations.

Set up all plates at the same time.

Track results over time

Compare outcomes, with respect to:

Same hits found?

Improvement (or not) in hits, labeled vs unlabeled?


Additional leads found (for labeled protein)?

1. Previously shown – Forsythe, Achari, & Pusey, Acta Cryst. D62, 339-346.

Diffraction Images – For Starters

IPPase, room temperature, home source, 10 min. exposure.

We have not found any effects on diffraction data from trace fluorescent labeling.



Crystals grown in and filled capillaries, mosaicity for both was equivalent

Ongoing Experiment Design/Plan

Test at least 30 different proteins – more if possible.

Soluble proteins – as available

Complexes – individual proteins and mixed

IMP's

Set up 6 plates for each protein – from the same stock solution/protein preparation.

3 labeled as per protocol

3 not labeled.

both at same final concentration.

use Nanodrop robot to set up plates

Use ONE screen – Hampton Research High Throughput 96 condition screen.

Follow plates for at least 8 weeks, scoring as described below.

Carry out optimization on all lead conditions (score = 4, to be described below).

Compare outcomes – labeled vs not.

Consistency Within Sets of Plates

A starting impetus for these experiments came from a review of six plates of one of the proteins (Tt141). A total of 38 hits were found (spheroids → 3D xtals), but not in all of the plates. These results represented at least two different protein preparations, with the plates set up at various times for each preparation.

Of the 38 total crystallization conditions found: (with still more testing, now up to ~50 conditions)

13 (34%) were found in all 6 plates – 1/3 of all conditions found in all 6 plates.

16 (42%) were found in 5 of the plates.

22 (58%) were found in 4 of the plates – 2/3 of the plates had just over 1/2 the conditions.


26 (68%) were found in 3 of the plates – 1/2 of the plates had just over 2/3 of the conditions.

30 (79%) were found in 2 of the plates.

NO ONE PLATE HAD ALL OF THE HITS FOUND! Two plates had 27 hits (71%), but did not have the same hits.

More Questions

Would I have found more conditions if I set up more plates using the same screen? What if you have 0 hits in a single plate w/a single screen? Is setting up more of the same plate approaching a definition for insanity?



Data Analysis

Scoring is on a (slightly) modified 9 point scale:

1 = clear solution

2 = phase separation

3 = precipitate

4 = precipitate with 'bright spots' – no discrete structure observed.

5 = non faceted structures (urchins, dendrites, spheroids) w/bright fluorescence – critical distinction vs. score of 4 is association with an observable structure.

6 = needles (pointed ends – flat ends = rods = 8 or 9)

7 = plate

8 = 3D crystals < 0.2 mm in size

9 = 3D crystals > 0.2 mm in size

Both labeled and unlabeled plates examined under white light for first scoring pass.


Labeled plates then examined under fluorescent illumination (primary source of score of 4) to adjust scores on basis of fluorescence image.

The well score is from the highest scoring structure in that well. Thus, if the well is mostly spheroids (score of 5), but one has developed facets, then the score is for a 3D crystal.


Unlabeled plate scores were 'edited' by reference to labeled plates. If similar structures appeared in both, but were not fluorescent, then the score was adjusted accordingly.

Scoring Examples


Score = 5




Score = 6




Score = 7



Score = 8



Score = 9



Score = 4

Examples to be shown...

Preliminary Results

By Conditions

17 proteins tested so far

Protein - activity

Protein - activity	-----SCORES OF-----				4 (leads) → Crystals	
	<u>5</u>	<u>6</u>	<u>7</u>	<u>8,9</u>	(Total/unique)	(Tested/found)
Tt55 – Translation Initiation Factor	26/21	12/14	3/1	3/3	31/14	14/4
Tt71 – Intracellular Protease	5/5	23/26	0/0	7/8	13/5	NT
Tt82 – HAD Superfamily Hydrolase	5/8	1/0	5/4	0/0	23/22	12/6
Tt81 – Haloacid Dehalogenase	42/54	9/5	0/0	18/16	11/0	NT
Tt80 – Sugar HAD Phosphatase	17/12	1/0	0/0	9/6	31/24	NT
Tt94 – RNA 3'-Terminal Phosphate Cycl.	9/8	11/11	0/0	8/8	21/14	NT
Tt102 – Endonuclease Methyltransferase	9/6	2/4	1/0	1/1	20/15	15/-IP
Tt106 – Nucleotide Kinase	5/3	0/0	0/0	0/0	16/15	14/4-IP ¹
Tt189 – Nucleoside Diphosphate Kinase	10/10	1/1	1/1	7/3	21/11	5/1
Tt186 – Alcohol Dehydrogenase	18/21	10/8	4/5	15/12	10/3	NT
Tt75 – Prolyl Endopeptidase	18/13	3/0	5/4	1/5	29/15	4/1
Tt46 – SS DNA-specif. Exonuclease	13/6	1/1	1/0	0/0	31/22	22/7
Tt141 – Inorganic Pyrophosphatase	21/23	13/14	4/6	20/12	11/3	1/1
Tt97 – Aspartate Racemase	24/2	3/2	0/0	6/5	14/8	NT
PCP – Pyrrolidone Carboxylate Peptidase	39/21	36/9	5/2	5/2	49/16	15/10-IP ¹
PCNA – Proliferating Cell Nuclear Antigen	22/22	10/6	3/0	14/11	21/8	NT
Phas – Kidney Bean Phaseolin	17/19	0/0	0/0	15/17	13/5	NT
	Totals 129/109				Ave. of <u>TESTED</u> 39 % (34/87)	

CRYSTALS FROM FLUORESCENCE


1) Capillary counter diffusion and IL-based optimization used.

Bright Spots in Precipitant For Lead Conditions

(Those outcomes with a score of 4)

Bright spots are often observed, which do not correspond to any identifiable discrete structures.


If intensity = structure, could they indicate proximity to crystallization conditions?



Canavalin – HSHT condition F1 (30%PEGMME 2000, 0.2M ammonium sulfate, 0.1M NaAcetate, pH 4.6).

Screening was first with respect to pH, then PEG MME & ammonium sulfate concentrations, and resulted in crystals as shown on the right (5% PEG MME 2000, 0.15M ammonium sulfate, 0.1M NaCacodylate, pH 6.5, 5 mg/ml protein). Scale bars are 200um.

Bright Spots in Precipitant For Lead Conditions





18% PEG 8K, 0.2M ZnAcetate,
0.1 M NaCacodylate, pH 6.5



Optimized to

10% PEG 8K, 0.01M ZnAcetate,
0.1M NaCacodylate, pH 6.5





β -Lactoglobulin B, left side shows results of HSHT screen condition D9.
This was optimized as shown on the right.

The Zn was necessary – crystals were not obtained without. Scale bars = 200 μ m.


Conclusion – Some precipitation conditions hid crystallization conditions.

More 'Bright Spot' Optimizations

Sometimes the spots are not all that bright




PCP, Hampton Research HT condition
C9. Optimization is condition C9 at
0.1 M C2[mim]-Cl, mixed 1:4
(v:v, screen:protein)



Scale bar = 200 um

Protein is Tt169, annotated as an exoribonuclease. The crystallizations conditions are an optimization of HSHT D6, the result of a 'bright spot' (below). The set-up is in a Microlytic counterdiffusion plate, with the imaged volume estimated as 15 nL.




**Benefits –
Disappearing salt crystals
the original reason for
using fluorescence**

**Native Canavalin,
labeled with CR**

HS II, #13



**Last night I saw upon the stair
A little man who wasn't there
He wasn't there again today
Oh, how I wish he'd go away.**



Another Example

(Joe thought he had me with this one!)

When the first diffraction attempt was made they centered the beam on the salt.




Finding Hidden Crystals

**Native Canavalin, N-terminal
amine labeled w/ the
fluorescent probe PyMPO**

HS 1, #14



**PyMPO was used in early
experiments as it had Ex and
Em characteristics similar to
GFP, for which the microscope
was set up.**





Finding Hidden Crystals

Native Canavalin, N-terminal
amine labeled w/ PyMPO

HS 1, #23





A Light Pipe Effect from Rod-Shaped Crystals



I have learned to closely examine precipitate showing pinpoints of light, which often indicate the presence of needles due to this light-pipe effect.

Quick Visual Scans are Facilitated Using Fluorescence Even at low resolution





Finding crystals can be a pain in the ass
Fluorescence is a method with real class
You really feel fine
Seeing them shine
And you can review the plates very fast!

Not recommended with excitation using UV light!

Future (now current) Directions

Multicolor Fluorescence for Complexes

Three different IPPase solutions individually labeled, each with a unique dye.
The solutions were mixed, then a screening plate set up.



From right to left-

Cascade Yellow, Ex @ 405 nm, Em 550-570 nm

Carboxyrhodamine, Ex @ 525 nm, Em 540-560 nm

Pacific Blue, Ex @ 405 nm, Em 450-460 nm.

Images acquired using original prototype system. This required changing the (single 5 mm) LED and the filter cube for each imaging run. The new system can do two colours and white light, or three colors without white light.

Multicolor Imaging



Protein is Tt141, HSHT condition E10

White Light




CR Source
& Filters

TR Source
& Filters




CR labeled
protein only



TR labeled
protein only



CR & TR labeled
protein




Probe selection is important! A higher wavelength (Ex and Em) red probe should be used instead of the TR. Probes chosen should be relatively immune from effects due to pH, solution components, etc.

Images acquired using X2 system.

Multicolor Fluorescence For Complexes – 1st Results

PCNA + Flap Endonuclease: The PCNA is labeled with Carboxyrhodamine (green) and the FEN with Texas Red (red). They were purified separately, mixed together at a 1:1 ratio, and a HSHT screening plate set up.



CAN IT BE USED WITH LCP?

Different labeling considerations when working with IMPs

Negatives

Fluorescent probes will usually (almost always) partition into the hydrophobic phase.

If not removed, detergent micelles can become a source of false positives.

Positives

Fluorescence yield is typically better in hydrophobic environments.

So – Fluorescence should be better for trace fluorescently labeled protein, but only if the free probe can be removed.

Solution So Far – Label the protein just after extraction and just prior to the first chromatography step. Use this step to do a detergent exchange, to remove unbound fluorescent probe. **First test did not work (very low yield of probe-protein product) and we are now rethinking this approach.**


Same considerations about pH and the presence of 1° and 2° amines applies (if labeling amines) when labeling.

Initial testing (and learning curve climbing) now underway with *E. coli* OmpF.

Future (Now Current) Directions An Automated Plate Imager

My Laboratory Prototype

~6 months assembly & programming



Has Evolved To This


(Thanks to Fitz-Thors Engineering!)






This prototype is still useful for testing out of new ideas, light sources, etc. However, It is limited in only being able to use a single wavelength. Routine data collection is now carried out using this system, 2 fluorescent and 1 white light colours (also possible to have 3 fluorescent colours, no white light). Eventual goal is 4 fluorescent colors & white light.

Future Directions – Automated Scoring



Intensity makes it very easy to identify objects of interest



B9a-5-22-50xf: Intensity Histogram



010504-CAN-HSB11a-32xf: Intensity Histogram



- Focus is important, although fluorescence is apparently more tolerant than white light images. Now Looking to apply autofocusing on images having sufficient contrast.
- False positives vs. false negatives – preferable to have 100% correct false negatives and 80% correct false positives than the reverse.
- Preliminary test goals were to divide the images into crystals, potential leads (bright spots that are not obvious crystals), precipitate, and clear.
- Crystals were successfully identified ~85% of the time, with missed crystals being primarily due to images out of focus, resulting in ‘bright spot’ lead conditions being classified as crystals and visa-versa.
- Future goals for image analysis software are to include morphological analysis, to facilitate differentiation of crystals from ‘bright spots’ and Include crystal morphology in the scoring.
- ‘Primitive’ scoring/outcome prioritization routine now built in to the Crystal X2 software.

Work in progress in collaboration with Dr. Ramazan Aygun, Dept. of Computer Sciences, UAH, as part of a Phase II STTR.

- Fluorescence can be a powerful tool in support of macromolecule crystallization, having many potential applications.
- It can be used to find lead crystallization conditions that would not be apparent using 'standard' plate-based methods.
- On a more 'routine' basis it can be used to rapidly find crystals or even potential lead conditions in plate screens.
- Intensity is a easier search parameter than straight lines, simplifying automated and/or visual plate-based screening analysis.
- Multi-color fluorescence can be used to verify the presence of multiple components in the crystallization of complexes .
- Probe derivatization levels of $\leq 0.5\%$ do not affect nucleation rates or crystal quality. – do not have to use modified protein for growth of crystals for analysis!! After all – you now have the Information you needed – the crystallization conditions.

Acknowledgements

iXpressGenes:

Dr. Joe Ng – crystallography

Michelle Morris – screening, crystal growth

Dr. Qunying Yuan – molecular biology

Jorge Barcena – protein purification

Univ. of Alabama, Huntsville

Dr. Ramazan Aygun – lead for image analysis software

Madhav Sigdel – student who does the work

Fitz-Thors Engineering, (recently moved to Birmingham, Alabama)

-Conversion of prototype system to the Crystal X2

Funding – NIH – Phase I & II STTRs (GM090453)

Any Questions?

