

Natrix™ 2 is a reagent kit designed to provide a rapid screening method for the crystallization of nucleic acids and nucleic acid-protein complexes. The screen is simple and practical for finding initial crystallization conditions as well as determining the solubility of nucleic acids in a wide range of precipitants and pH.

The kit is designed to provide a biased sparse matrix of trial conditions selected from known and published crystallization conditions. The reagent parameter variables are pH, buffer, salt, additive, and precipitant. Six different pH's: 5.5, 6.0, 6.5, 7.0, 7.5, and 8.5 are utilized with Sodium cacodylate, MOPS, HEPES sodium, PIPES, and TRIS hydrochloride as the buffers. The four categories of precipitating agents utilized are volatile agents, non-volatile agents, salts, and a combination of these three. Refer to the enclosed Natrix 2 reagent formulation for additional information.

Sample Preparation

The sample should be as pure as is practically possible (> 95%) and free of amorphous and particulate material. Remove amorphous material by centrifugation (or micro-filtration when appropriate) prior to use.

Recommended stock concentration of the nucleic acid is 0.5 to 1.0 mM or 5 to 10 mg/ml depending upon the solubility and size of the sample. The nucleic acid should be solubilized in a water based system which promotes the stability and monodispersity of the nucleic acid. If a buffer is utilized for nucleic acid preparation, a concentration of 5 to 10 mM is recommended in order to allow the buffers in Natrix 2 to alter the pH of the sample drop.

One may wish to include a polyamine such as spermine or spermidine at a concentration of 0.5 to 1.5 mM. The polyamine need not be added to the reservoir. Finally, when sample annealing is desired, the sample should be preheated to 50°C to 95°C for 10 minutes then cooled slowly to 25°C (room temperature), in the presence of the samples buffer and 5-20 mM Magnesium chloride to produce native molecules.¹ After cooling, centrifuge and micro-filtrate the sample.


Performing The Screen

The following procedure describes the use of Natrix 2 with the Hanging Drop Vapor Diffusion method. Natrix 2 is also compatible with the Sitting Drop, Sandwich Drop, Microbatch, Microdialysis, and Free Interface Diffusion methods. A complete description of the Hanging, Sitting, Sandwich Drop, Dialysis and other crystallization methods are available from the Hampton Research Crystal Growth 101 Library.

1. Prepare a VDX Plate (HR3-140) for Hanging Drop Vapor Diffusion by applying a thin bead of cover slide sealant to the upper edge of each of the 24 reservoirs. One may also use a Greased VDX Plate (HR3-170). Forty eight reservoirs are to be prepared for a complete Natrix 2. See Figure 1.

Figure 1

Cross section of a reservoir in the VDX plate.



2. Using a clean pipet tip, pipet 1 ml of Natrix 2 reagent 1 into reservoir A1. Discard the pipet tip, add a new pipet tip and pipet 1 ml of Natrix 2 reagent 2 into reservoir A2. Repeat the procedure for the remaining 46 Natrix 2 reagents using a clean pipet tip for each reagent so as to avoid reagent contamination and carry over.

Figure 2

Crystallization Droplet (2 µl Sample / 2 µl Reagent)

Siliconized Coverslip


3. Pipet 2 µl of the sample to the center of a clean, siliconized 22 mm diameter circle or square cover slide. See Figure 2.

4. Pipet 2 µl of Natrix 2 reagent 1 from reservoir A1 into the sample droplet and mix by aspirating and dispensing the droplet several times, keeping the tip in the drop during mixing to avoid foaming. See Figure 2.

5. Working quickly to minimize evaporation, invert the cover slide and droplet over reservoir A1 and seal the cover slide onto the edge of the reservoir. See Figure 3.

Figure 3

Inverted siliconized coverslip placed over the reservoir.




6. Repeat operations 3 through 5 for the remaining 47 Natrix 2 reagents.

7. If the quantity of sample permits, perform Natrix 2 in duplicate and incubate one set of plates at 4°C and the second set at room temperature. Incubate and store the crystallization plates in a stable temperature environment free of vibration.

Examine The Drop

Carefully examine the drops under a stereo microscope (10 to 100x magnification) immediately after setting up the screen. Record all observations and be particularly careful to scan the focal plane for small crystals. Observe the drops once each day for the first week, then once a week thereafter. Records should indicate whether the drop is clear, contains precipitate, and/or crystals. It is helpful to describe the drop contents using descriptive terms.

Figure 4
Typical observations in a crystallization experiment



Adding magnitude is also helpful. Example: 4+ yellow/brown fine precipitate, 2+ small bipyramid crystals, clear drop, 3+ needle shaped crystals in 1+ white precipitate. One may also employ a standard numerical scoring scheme (Clear = 0, Precipitate = 1, Crystal = 10, etc). Figure 4 (on page 2) shows typical examples of what one might observe in a crystallization experiment.

Interpreting Matrix 2

Clear drops indicate that either the relative supersaturation of the sample and reagent is too low or the drop has not yet completed equilibration. If the drop remains clear after 3 to 4 weeks consider repeating the Matrix 2 condition and doubling the sample concentration. If more than 33 of the 48 Matrix 2 drops are clear consider doubling the sample concentration and repeating the entire screen.

Drops containing precipitate indicate that either the relative supersaturation of the sample and reagent is too high, the sample has denatured, or the sample is heterogeneous. To reduce the relative supersaturation, dilute the sample twofold and repeat the Matrix 2 condition. If more than 33 of the 48 Matrix 2 drops contain precipitate and no crystals are present, consider diluting the sample concentration in half and repeating the entire screen. If sample denaturation is suspect, take measures to stabilize the sample (add reducing agent, ligands, glycerol, salt, or other stabilizing agents). If the sample is impure, aggregated, or heterogeneous take measures to pursue homogeneity. It is possible to obtain crystals from precipitate so do not discard nor ignore a drop containing precipitate. If possible, examine drops containing precipitate under polarizing optics to differentiate precipitate from microcrystalline material.

If the drop contains a macromolecular crystal the relative supersaturation of the sample and reagent is good. The next step is to optimize the preliminary conditions (pH, salt type, salt concentration, precipitant type, precipitant concentration, sample concentration, temperature, additives, and other crystallization variables) which produced the crystal in order to improve crystal size and quality.

Compare the observations between the 4°C and room temperature incubation to determine the effect of temperature on sample solubility. Different results in the same drops at different temperatures indicate that sample solubility is temperature dependent and that one should include temperature as a

variable in subsequent screens and optimization experiments.

Retain and observe plates until the drops are dried out. Crystal growth can occur within 15 minutes or one year.

Matrix 2 Formulation

Matrix 2 reagents are formulated using the highest purity chemicals, ultrapure water (18.2 Megohm-cm, 5 ppb TOC) and are sterile filtered using 0.22 micron filters into sterile containers (no preservatives added).

Matrix 2 reagents are readily reproduced using Hampton Research Optimize™ stock solutions of salts, polymers and buffers. Optimize stock reagents make reproducing Matrix 2 reagents fast, convenient and easy. Dilutions can be performed directly into the crystallization plate using Optimize stock reagents.

Matrix 2 reagents containing buffers are formulated by creating a 1.0 M stock buffer, titrated to the desired pH using Hydrochloric acid or Sodium hydroxide. The buffer is then diluted with the other reagent components and water. No further pH adjustment is required.

Matrix 2 reagents are stable at room temperature and are best if used within 12 months of receipt. To enhance reagent stability it is recommended that Matrix 2 be stored at 4°C or -20°C. Avoid ultraviolet light to preserve reagent stability.

If the sample contains phosphate, borate, or carbonate buffers it is possible to obtain inorganic crystals (false positives) when using Matrix 2 reagents containing divalent cations such as magnesium, calcium, or zinc. To avoid false positives use phosphate, borate, or carbonate buffers at concentrations of 10 mM or less or exchange the phosphate, borate, or carbonate buffer with a more soluble buffer that does not complex with divalent cations.

References and Readings

1. Crystallization of ribozymes and small RNA motifs by a sparse matrix approach. Doudna, J.A. et al. PNAS, USA, 90, 7829-7833, 1993.
2. Doudna, J.A., Grosshans, C., Godding, A., and Kundrot, C.E., Proc. Natl. Acad. Sci. USA, Aug 1993, Vol. 90, 7829-7833.

3. McPherson, A., Current Approaches to Macromolecular Crystallization, Eur. J. Biochem., 189, 1-23, 1990.
4. Protein and Nucleic Acid Crystallization, Methods: A Companion to Methods in Enzymology, Academic Press, Vol 1, No. 1, Aug 1990.
5. Crystallization of Nucleic Acids and Proteins - The Practical Approach Series, Edited by A. Ducruix and R. Giege, Oxford Univ. Press 1992.
6. Cudney, R., et al, Acta Cryst, D50, 414-423, 1994.
7. Schultz, S.C., Shields, G.C., and Steitz, T.A., J. of Mol. Bio., 213, 159-166, 1990.
8. Baeyens, K.J., Jancarik, J., and Holbrook, S.R., Acta Cryst, D50, 764-767, 1994.
9. Perona, J.J., Methods: A Companion to Methods in Enzymology, Vol. 1, 75-82, 1990.
10. Oubridge, C., Ito, N., Teo, C.-H., Fearnley, I., and Nagai, K., J. Mol. Bio., 249, 408-412, 1995.
11. Price, S.R., Ito, N., Oubridge, C., Avis, J.M., and Nagai, K., J. Mol. Bio., 249, 398-408.
12. Wincott, F., DiRenzo, A., Shaffer, C., Scaringe, S., and Usman, N., Nucleic Acid Res., Vol. 23, 2677-2687, 1995.

Technical Support

Inquiries regarding Matrix 2 reagent formulation, interpretation of screen results, optimization strategies and general inquiries regarding crystallization are welcome. Please e-mail, fax, or telephone your request to Hampton Research. Fax and e-mail Technical Support are available 24 hours a day. Telephone technical support is available 8:00 a.m. to 4:30 p.m. USA Pacific Standard Time

Hampton Research
34 Journey
Aliso Viejo, CA 92656-3317 U.S.A.
Tel: (949) 425-1321 • Fax: (949) 425-1611
Technical Support e-mail: tech@hrmail.com
Website: www.hamptonresearch.com

Tube #	Salt	Tube #	Buffer ◇	Tube #	Precipitant	Tube #	Additive
1.	0.04 M Lithium chloride, 0.02 M Magnesium chloride hexahydrate	1.	0.04 M Sodium cacodylate trihydrate pH 5.5	1.	30% v/v (+/-)-2-Methyl-2,4-pentanediol	1.	0.002 M Hexamine cobalt(III) chloride
2.	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate	2.	0.04 M Sodium cacodylate trihydrate pH 5.5	2.	35% v/v (+/-)-2-Methyl-2,4-pentanediol	2.	0.002 M Hexamine cobalt(III) chloride
3.	0.012 M Sodium chloride, 0.08 M Potassium chloride	3.	0.04 M Sodium cacodylate trihydrate pH 5.5	3.	45% v/v (+/-)-2-Methyl-2,4-pentanediol	3.	0.002 M Hexamine cobalt(III) chloride
4.	0.02 M Magnesium chloride hexahydrate	4.	0.04 M Sodium cacodylate trihydrate pH 5.5	4.	40% v/v (+/-)-2-Methyl-2,4-pentanediol	4.	0.002 M Hexamine cobalt(III) chloride
5.	0.002 M Calcium chloride dihydrate	5.	0.05 M Sodium cacodylate trihydrate pH 6.0	5.	1.8 M Ammonium sulfate	5.	0.0005 M Spermine
6.		6.	0.05 M Sodium cacodylate trihydrate pH 6.0	6.	35% v/v Tacsimate™ pH 6.0	6.	0.001 M Spermine
7.	0.1 M Sodium chloride	7.	0.05 M Sodium cacodylate trihydrate pH 6.0	7.	10% w/v Polyethylene glycol 4,000	7.	0.0005 M Spermine
8.	0.05 M Potassium chloride	8.	0.05 M Sodium cacodylate trihydrate pH 6.0	8.	10% w/v Polyethylene glycol 8,000	8.	0.0005 M Spermine, 0.0005 M L-Argininamide dihydrochloride
9.	0.1 M Potassium chloride	9.	0.05 M Sodium cacodylate trihydrate pH 6.0	9.	16% w/v Polyethylene glycol 1,000	9.	0.0005 M Spermine
10.	0.005 M Magnesium chloride hexahydrate, 0.002 M Calcium chloride dihydrate	10.	0.05 M Sodium cacodylate trihydrate pH 6.0	10.	15% v/v 2-Propanol	10.	0.001 M Spermine
11.	0.075 M Sodium chloride, 0.002 M Calcium chloride dihydrate	11.	0.05 M Sodium cacodylate trihydrate pH 6.0	11.	30% w/v 1,6-Hexanediol	11.	0.0005 M Spermine
12.	0.02 M Magnesium sulfate hydrate, 0.002 M Cobalt(II) chloride hexahydrate	12.	0.05 M Sodium cacodylate trihydrate pH 6.0	12.	25% v/v (+/-)-2-Methyl-2,4-pentanediol	12.	0.0005 M Spermine
13.		13.	0.05 M Sodium cacodylate trihydrate pH 6.0	13.	30% v/v (+/-)-2-Methyl-2,4-pentanediol	13.	
14.	0.08 M Sodium chloride, 0.012 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	14.	0.04 M Sodium cacodylate trihydrate pH 6.0	14.	30% v/v (+/-)-2-Methyl-2,4-pentanediol	14.	0.012 M Spermine tetrahydrochloride
15.	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate	15.	0.04 M Sodium cacodylate trihydrate pH 6.0	15.	35% v/v (+/-)-2-Methyl-2,4-pentanediol	15.	0.012 M Spermine tetrahydrochloride
16.	0.08 M Strontium chloride hexahydrate	16.	0.04 M Sodium cacodylate trihydrate pH 6.0	16.	35% v/v (+/-)-2-Methyl-2,4-pentanediol	16.	0.012 M Spermine tetrahydrochloride
17.	0.08 M Potassium chloride, 0.02 M Barium chloride dihydrate	17.	0.04 M Sodium cacodylate trihydrate pH 6.0	17.	40% v/v (+/-)-2-Methyl-2,4-pentanediol	17.	0.012 M Spermine tetrahydrochloride
18.	0.08 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	18.	0.04 M Sodium cacodylate trihydrate pH 6.0	18.	45% v/v (+/-)-2-Methyl-2,4-pentanediol	18.	0.012 M Spermine tetrahydrochloride
19.	0.08 M Sodium chloride	19.	0.04 M Sodium cacodylate trihydrate pH 6.0	19.	45% v/v (+/-)-2-Methyl-2,4-pentanediol	19.	0.012 M Spermine tetrahydrochloride
20.	0.08 M Sodium chloride, 0.02 M Barium chloride dihydrate	20.	0.04 M Sodium cacodylate trihydrate pH 6.0	20.	45% v/v (+/-)-2-Methyl-2,4-pentanediol	20.	0.012 M Spermine tetrahydrochloride
21.	0.012 M Sodium chloride, 0.08 M Potassium chloride	21.	0.04 M Sodium cacodylate trihydrate pH 6.0	21.	50% v/v (+/-)-2-Methyl-2,4-pentanediol	21.	0.012 M Spermine tetrahydrochloride
22.	0.08 M Potassium chloride	22.	0.04 M Sodium cacodylate trihydrate pH 6.0	22.	55% v/v (+/-)-2-Methyl-2,4-pentanediol	22.	0.012 M Spermine tetrahydrochloride
23.	0.018 M Magnesium chloride hexahydrate	23.	0.05 M Sodium cacodylate trihydrate pH 6.5	23.	10% v/v 2-Propanol	23.	0.003 M Spermine
24.	0.02 M Magnesium chloride hexahydrate	24.	0.05 M MOPS pH 7.0	24.	2.0 M Ammonium sulfate	24.	0.0005 M Spermine

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

Natrix 2 contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.

Tube #	Salt	Tube #	Buffer ◊	Tube #	Precipitant	Tube #	Additive
25.		25.	0.05 M HEPES sodium pH 7.0	25.	40% v/v Tacsimate™ pH 7.0	25.	0.002 M Spermine, 0.002 M Hexamine cobalt(III) chloride
26.	0.02 M Magnesium chloride hexahydrate	26.	0.05 M MOPS pH 7.0	26.	55% v/v Tacsimate™ pH 7.0	26.	0.002 M Hexamine cobalt(III) chloride
27.	0.02 M Magnesium chloride hexahydrate	27.	0.05 M Sodium cacodylate trihydrate pH 7.0	27.	15% v/v 2-Propanol	27.	0.001 M Hexamine cobalt(III) chloride, 0.001 M Spermine
28.	0.005 M Magnesium chloride hexahydrate	28.	0.05 M MOPS pH 7.0	28.	25% v/v 1,4-Dioxane	28.	0.001 M Spermine
29.	0.01 M Magnesium chloride hexahydrate, 0.002 M Barium chloride dihydrate	29.	0.05 M MOPS pH 7.0	29.	30% v/v 1,4-Dioxane	29.	
30.	0.001 M Magnesium chloride hexahydrate, 0.002 M Calcium chloride dihydrate	30.	0.05 M MOPS pH 7.0	30.	15% v/v (+/-)-2-Methyl-2,4-pentanediol	30.	
31.	0.08 M Strontium chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate	31.	0.04 M Sodium cacodylate trihydrate pH 7.0	31.	20% v/v (+/-)-2-Methyl-2,4-pentanediol	31.	0.012 M Spermine tetrahydrochloride
32.	0.08 M Sodium chloride	32.	0.04 M Sodium cacodylate trihydrate pH 7.0	32.	30% v/v (+/-)-2-Methyl-2,4-pentanediol	32.	0.012 M Spermine tetrahydrochloride
33.	0.04 M Lithium chloride, 0.08 M Strontium chloride hexahydrate	33.	0.04 M Sodium cacodylate trihydrate pH 7.0	33.	30% v/v (+/-)-2-Methyl-2,4-pentanediol	33.	0.012 M Spermine tetrahydrochloride
34.	0.04 M Lithium chloride, 0.08 M Strontium chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate	34.	0.04 M Sodium cacodylate trihydrate pH 7.0	34.	30% v/v (+/-)-2-Methyl-2,4-pentanediol	34.	0.012 M Spermine tetrahydrochloride
35.	0.08 M Sodium chloride, 0.012 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	35.	0.04 M Sodium cacodylate trihydrate pH 7.0	35.	35% v/v (+/-)-2-Methyl-2,4-pentanediol	35.	0.012 M Spermine tetrahydrochloride
36.	0.012 M Sodium chloride, 0.08 M Potassium chloride	36.	0.04 M Sodium cacodylate trihydrate pH 7.0	36.	40% v/v (+/-)-2-Methyl-2,4-pentanediol	36.	0.012 M Spermine tetrahydrochloride
37.	0.08 M Sodium chloride, 0.02 M Barium chloride dihydrate	37.	0.04 M Sodium cacodylate trihydrate pH 7.0	37.	40% v/v (+/-)-2-Methyl-2,4-pentanediol	37.	0.012 M Spermine tetrahydrochloride
38.	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate	38.	0.04 M Sodium cacodylate trihydrate pH 7.0	38.	40% v/v (+/-)-2-Methyl-2,4-pentanediol	38.	0.012 M Spermine tetrahydrochloride
39.	0.08 M Potassium chloride, 0.02 M Barium chloride dihydrate	39.	0.04 M Sodium cacodylate trihydrate pH 7.0	39.	40% v/v (+/-)-2-Methyl-2,4-pentanediol	39.	0.012 M Spermine tetrahydrochloride
40.	0.08 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	40.	0.04 M Sodium cacodylate trihydrate pH 7.0	40.	50% v/v (+/-)-2-Methyl-2,4-pentanediol	40.	0.012 M Spermine tetrahydrochloride
41.	0.08 M Potassium chloride	41.	0.04 M Sodium cacodylate trihydrate pH 7.0	41.	60% v/v (+/-)-2-Methyl-2,4-pentanediol	41.	0.012 M Spermine tetrahydrochloride
42.	0.02 M Magnesium chloride hexahydrate, 0.002 M Cobalt(II) chloride hexahydrate	42.	0.05 M HEPES sodium pH 7.5	42.	2.0 M Ammonium sulfate	42.	0.001 M Spermine
43.	0.02 M Magnesium chloride hexahydrate	43.	0.05 M PIPES pH 7.5	43.	4% w/v Polyethylene glycol 8,000	43.	0.001 M Spermine
44.	0.015 M Magnesium chloride hexahydrate, 0.002 M Barium chloride dihydrate	44.	0.05 M PIPES pH 7.5	44.	7% v/v 2-Propanol	44.	0.0005 M Spermine
45.	0.02 M Magnesium chloride hexahydrate	45.	0.05 M PIPES pH 7.5	45.	10% w/v 1,6-Hexanediol	45.	0.001 M Spermine
46.	0.01 M Magnesium chloride hexahydrate	46.	0.05 M HEPES sodium pH 7.5	46.	15% v/v (+/-)-2-Methyl-2,4-pentanediol	46.	0.0015 M Spermine
47.	0.2 M Calcium chloride dihydrate	47.	0.05 M HEPES sodium pH 7.5	47.	28% v/v Polyethylene glycol 400	47.	0.002 M Spermine
48.	0.002 M Copper(II) chloride dihydrate	48.	0.05 M TRIS hydrochloride pH 8.5	48.	1.8 M Lithium sulfate monohydrate	48.	0.0005 M Spermine

◊ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

34 Journey
 Aliso Viejo, CA 92656-3317 U.S.A.
 Tel: (949) 425-1321 • Fax: (949) 425-1611
 E-mail: tech@hrmail.com
 Website: www.hamptonresearch.com

Natrix 2 contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.



Solutions for Crystal Growth

© 1991 - 2021 Hampton Research Corp. all rights reserved
 Printed in the United States of America. This guide or
 parts thereof may not be reproduced in any form without
 the written permission of the publishers.

Sample: _____ Sample Concentration: _____
 Sample Buffer: _____ Date: _____
 Reservoir Volume: _____ Temperature: _____
 Drop Volume: Total _____ μ l Sample _____ μ l Reservoir _____ μ l Additive _____ μ l

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals
- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2 mm)
- 9 Single Crystals (3D Growth > 0.2 mm)

Matrx™ 2 - HR2-117 Scoring Sheet

Date: _____ Date: _____

1.	0.04 M Lithium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 5.5, 30% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride		
2.	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 5.5, 35% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride		
3.	0.012 M Sodium chloride, 0.08 M Potassium chloride, 0.04 M Sodium cacodylate trihydrate pH 5.5, 45% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride		
4.	0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 5.5, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride		
5.	0.002 M Calcium chloride dihydrate, 0.05 M Sodium cacodylate trihydrate pH 6.0, 1.8 M Ammonium sulfate, 0.0005 M Spermine		
6.	0.05 M Sodium cacodylate trihydrate pH 6.0, 35% v/v Tacsimate pH 6.0, 0.001 M Spermine		
7.	0.1 M Sodium chloride, 0.05 M Sodium cacodylate trihydrate pH 6.0, 10% w/v Polyethylene glycol 4,000, 0.0005 M Spermine		
8.	0.05 M Potassium chloride, 0.05 M Sodium cacodylate trihydrate pH 6.0, 10% w/v Polyethylene glycol 8,000, 0.0005 M Spermine, 0.0005 M L-Argininamide dihydrochloride		
9.	0.1 M Potassium chloride, 0.05 M Sodium cacodylate trihydrate pH 6.0, 16% w/v Polyethylene glycol 1,000, 0.0005 M Spermine		
10.	0.005 M Magnesium chloride hexahydrate, 0.002 M Calcium chloride dihydrate, 0.05 M Sodium cacodylate trihydrate pH 6.0, 15% v/v 2-Propanol, 0.001 M Spermine		
11.	0.075 M Sodium chloride, 0.002 M Calcium chloride dihydrate, 0.05 M Sodium cacodylate trihydrate pH 6.0, 30% w/v 1,6-Hexanediol, 0.0005 M Spermine		
12.	0.02 M Magnesium sulfate hydrate, 0.002 M Cobalt(II) chloride hexahydrate, 0.05 M Sodium cacodylate trihydrate pH 6.0, 25% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.0005 M Spermine		
13.	0.05 M Sodium cacodylate trihydrate pH 6.0, 30% v/v (+/-)-2-Methyl-2,4-pentanediol		
14.	0.08 M Sodium chloride, 0.012 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 6.0, 30% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
15.	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 6.0, 35% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
16.	0.08 M Strontium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 6.0, 35% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
17.	0.08 M Potassium chloride, 0.02 M Barium chloride dihydrate, 0.04 M Sodium cacodylate trihydrate pH 6.0, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
18.	0.08 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 6.0, 45% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
19.	0.08 M Sodium chloride, 0.04 M Sodium cacodylate trihydrate pH 6.0, 45% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
20.	0.08 M Sodium chloride, 0.02 M Barium chloride dihydrate, 0.04 M Sodium cacodylate trihydrate pH 6.0, 45% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
21.	0.012 M Sodium chloride, 0.08 M Potassium chloride, 0.04 M Sodium cacodylate trihydrate pH 6.0, 50% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
22.	0.08 M Potassium chloride, 0.04 M Sodium cacodylate trihydrate pH 6.0, 55% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
23.	0.018 M Magnesium chloride hexahydrate, 0.05 M Sodium cacodylate trihydrate pH 6.5, 10% v/v 2-Propanol, 0.003 M Spermine		
24.	0.02 M Magnesium chloride hexahydrate, 0.05 M MOPS pH 7.0, 2.0 M Ammonium sulfate, 0.0005 M Spermine		
25.	0.05 M HEPES sodium pH 7.0, 40% v/v Tacsimate pH 7.0, 0.002 M Spermine, 0.002 M Hexamine cobalt(III) chloride		
26.	0.02 M Magnesium chloride hexahydrate, 0.05 M MOPS pH 7.0, 55% v/v Tacsimate pH 7.0, 0.002 M Hexamine cobalt(III) chloride		
27.	0.02 M Magnesium chloride hexahydrate, 0.05 M Sodium cacodylate trihydrate pH 7.0, 15% v/v 2-Propanol, 0.001 M Hexamine cobalt(III) chloride, 0.001 M Spermine		
28.	0.005 M Magnesium chloride hexahydrate, 0.05 M MOPS pH 7.0, 25% v/v 1,4-Dioxane, 0.001 M Spermine		
29.	0.01 M Magnesium chloride hexahydrate, 0.002 M Barium chloride dihydrate, 0.05 M MOPS pH 7.0, 30% v/v 1,4-Dioxane		
30.	0.001 M Magnesium chloride hexahydrate, 0.002 M Calcium chloride dihydrate, 0.05 M MOPS pH 7.0, 15% v/v (+/-)-2-Methyl-2,4-pentanediol		
31.	0.08 M Strontium chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 20% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
32.	0.08 M Sodium chloride, 0.04 M Sodium cacodylate trihydrate pH 7.0, 30% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
33.	0.04 M Lithium chloride, 0.08 M Strontium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 30% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
34.	0.04 M Lithium chloride, 0.08 M Strontium chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 30% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
35.	0.08 M Sodium chloride, 0.012 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 35% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
36.	0.012 M Sodium chloride, 0.08 M Potassium chloride, 0.04 M Sodium cacodylate trihydrate pH 7.0, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
37.	0.08 M Sodium chloride, 0.02 M Barium chloride dihydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
38.	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
39.	0.08 M Potassium chloride, 0.02 M Barium chloride dihydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
40.	0.08 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 7.0, 50% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
41.	0.08 M Potassium chloride, 0.04 M Sodium cacodylate trihydrate pH 7.0, 60% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.012 M Spermine tetrahydrochloride		
42.	0.02 M Magnesium chloride hexahydrate, 0.002 M Cobalt(II) chloride hexahydrate, 0.05 M HEPES sodium pH 7.5, 2.0 M Ammonium sulfate, 0.001 M Spermine		
43.	0.02 M Magnesium chloride hexahydrate, 0.05 M PIPES pH 7.5, 4% w/v Polyethylene glycol 8,000, 0.001 M Spermine		
44.	0.015 M Magnesium chloride hexahydrate, 0.002 M Barium chloride dihydrate, 0.05 M PIPES pH 7.5, 7% v/v 2-Propanol, 0.0005 M Spermine		
45.	0.02 M Magnesium chloride hexahydrate, 0.05 M PIPES pH 7.5, 10% w/v 1,6-Hexanediol, 0.001 M Spermine		
46.	0.01 M Magnesium chloride hexahydrate, 0.05 M HEPES sodium pH 7.5, 15% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.0015 M Spermine		
47.	0.2 M Calcium chloride dihydrate, 0.05 M HEPES sodium pH 7.5, 28% v/v Polyethylene glycol 400, 0.002 M Spermine		
48.	0.002 M Copper(II) chloride dihydrate, 0.05 M TRIS hydrochloride pH 8.5, 1.8 M Lithium sulfate monohydrate, 0.0005 M Spermine		