

Applications

Preliminary screen for the crystallization of membrane proteins as well as soluble biological macromolecules.

Features

- Single 96 Deep Well block format
- Compatible with robotics and multi-channel pipets
- Primary screen variables : salt, pH, and precipitant (salts, polymers, volatile organics, and non-volatile organics)
- Samples pH 4.6 – 8.5
- Membrane protein sparse matrix screen reagent conditions based on MembFac™ and Crystal Screen Lite™
- Preformulated, ready to screen
- Formulated for use with detergents

General Description

MembFac HT™ is a highly effective sparse matrix specifically designed as a preliminary screen for the crystallization of membrane proteins as well as soluble biological macromolecules. MembFac HT is designed as a 96 reagent crystallization screen that combines the strategies of MembFac and Crystal Screen Lite into a highly effective and efficient format. This kit allows one to evaluate a large variety of potential crystallization conditions with the 96 unique reagents.

MembFac HT is supplied in a sterile, polypropylene 96 Deep Well block, each reservoir containing 1 ml of sterile filtered reagent. The block is compatible with robotic and multi-channel pipet liquid handling systems and is heat sealed using a special polypropylene backed film. Each MembFac HT kit is supplied with an adhesive sealing film which can be used to seal the block after removing the heat seal. Additional adhesive sealing films can be obtained from Hampton Research or laboratory supply companies which offer high throughput plates and seals.

Within the 96 Deep Well block, rows A through D feature the 48 reagents of MembFac (HR2-114). These reagent variables consist of pH, buffer material salt, salt, and precipitant. Five different pH's: 4.6, 5.6, 6.5, 7.5, and 8.5 are utilized with Sodium acetate trihydrate, Sodium citrate tribasic dihydrate, ADA, HEPES sodium, and TRIS hydrochloride as the reagent buffers. The four categories of precipitating reagents consist of: volatile agents, non-volatile agents, salts, and a combination of these three. These reagents were selected specifically for use with detergents.

Rows E through H feature the first 48 reagents of Crystal Screen Lite (HR2-128). Based upon the original Jancarik and Kim screen also known as Crystal Screen and designed to provide a rapid screening method for the crystallization of biological macromolecules, this kit is also effective in determining the solubility of these macromolecules in a wide range of precipitants and

pH.³ The primary screens are salt, pH, and precipitant (salts, polymers, volatile organics, and non-volatile organics).² Crystal Screen Lite differs from the original Crystal Screen kit such that the primary precipitant reagents are one-half the concentration of that used in the original screen. The secondary salts, ion, and buffers remain at the original Crystal Screen concentration. Reducing the primary concentration of the primary precipitant results in a screen which is “more gentle” on the sample and typically produces much less precipitate conditions than the original Crystal Screen. Results comparing the Crystal Screen Lite formulation versus simply diluting the Crystal Screen formulation two-fold demonstrated more crystals using the Crystal Screen Lite protocol than the two-fold diluted Crystal Screen illustrating the importance of retaining the original salt, ion, and buffer concentration in Crystal Screen.⁵ Results comparing simply diluting the sample versus using Crystal Screen Lite also demonstrated more crystals when using Crystal Screen Lite than when simply diluting the sample. Crystal Screen Lite should be used with samples which demonstrate limited solubility in traditional crystallization reagents.^{1,4}

Refer to the enclosed MembFac HT reagent formulation for additional information on all 96 reagents.

Sample Preparation

The membrane protein of interest is isolated in the detergent which gives the highest stability/activity ratio. The final protein concentration should be 10 to 20 mg/ml and the detergent concentration should only be slightly above the CMC.

The sample should be as pure as is practically possible (>95%) and free of amorphous and particulate material. Remove amorphous material by centrifugation prior to use.

For additional sample preparation recommendation see Crystal Growth 101 - Preliminary Sample Preparation bulletin from Hampton Research.

Preparing the Deep Well Block for Use


It is recommended the Deep Well block be centrifuged before removing the sealing film. Centrifugation at 500 rpm for five minutes will remove stray reagent from the sealing film. Removing the reagent from the film prevents stray reagent droplets from falling into neighboring wells during film removal. After centrifugation the film can be removed by grasping a corner of the film and gently peeling the film from the plate. Alternatively, the film can be left intact and then pierced for reagent access.

Performing the Screen

Manual Method - Sitting Drop Vapor Diffusion


1. Using a 96 well sitting drop vapor diffusion plate, pipet the recommended

Figure 6
Typical observations in a crystallization experiment




volume (typically 50 to 100 microliters) of crystallization reagent from the Deep Well block into the reservoirs of the crystallization plate. The Deep Well block is compatible with 8 and 12 channel pipets as well as many automated liquid handling systems. Use clean pipet tips for each reagent set transfer and change pipet tips when changing reagents. For an 8 channel pipet, transfer reagents A1-H1 to reservoirs A1-H1 of the crystallization plate. Repeat this procedure for reagent columns B through H. Change pipet tips when moving between reagent columns. For a 12 channel pipet, transfer reagents A1-A12 to reservoirs A1-A12 of the crystallization plate. Repeat this procedure for reagent rows 1 through 12. See Figure 1. Time and pipet tips can be conserved by batch pipetting multiple plates with the same (row or column) of reagent before changing reagent and pipet tips.

Figure 1



2. Using clean pipet tips, pipet 0.05 to 2 microliters of crystallization reagent from the crystallization plate reservoir to the sitting drop well. Some 96 well crystallization plates allow this procedure to be performed using a multichannel pipet where other plates require the use of a single channel pipet. Change the pipet tip between reagents. See Figure 2.

Figure 2



3. Using a clean pipet tip, pipet 0.05 to 2 microliters of sample to the reagent drop in the sitting drop well. One may choose to simply dispense the sample with no mixing or dispense with mixing by gently aspirating and dispensing the sample several times, keeping the tip in the drop during mixing to avoid foaming. Work carefully but quickly to minimize evaporation from the crystallization plate. See Figure 2 above.

4. Seal the crystallization plate as per the manufacturer's recommendation. Most 96 well crystallization plates are sealed


using a clear sealing tape or film. View and score the experiment as desired. See Hampton Research technical bulletin Crystal Growth 101 - Viewing Crystallization Experiments for additional information on viewing drops.

5. Seal the remaining reagent in the Deep Well block using either clear sealing tape, film, or cap mat.

Manual Method – Microbatch 96 well format


1. Using a 96 well clear polystyrene microplate (U-bottom recommended for best drop centering, flat-bottom recommended for best optics) pipet 50 to 150 microliters of microbatch compatible oil into each of the 96 reservoirs. This can be accomplished using an 8 or 12 channel pipet and pipetting the oil from a reagent basin. See Figure 3.

Figure 3




2. Once the plate is oiled, use an 8 or 12 channel pipet to aspirate reagent from the Deep Well block and dispense the reagent under the oil in the Microbatch plate. Change tips when changing reagent to prevent cross reagent contamination. To save time and pipet tips, set multiple plates at one time. See Figure 4.

Figure 4



3. Using a single channel pipet, aspirate the sample and dispense the sample under oil in the Microbatch plate. It is not necessary to dispense the sample drop into the reagent drop or mix the drops. See Figure 5.

Figure 5



4. After all reagent and sample drops have been dispensed to the Microbatch plate, place the loose fitting clear cover on the Microbatch plate and centrifuge the plate for 10 minutes at 500 rpm. Centrifugation will cause the drops to coalesce into a single drop.

Note: If the drops appear flat or is fragmented into multiple drops, the centrifugation speed is too high and the centrifugation time is too long - adjust to obtain a spherical single drop in the center of the well.

5. Store the plates with the loose fitting clear polystyrene cover and observe for crystals. See Hampton Research technical bulletin Crystal Growth 101 - Viewing Crystallization Experiments for additional information on viewing drops.

MembFac HT Deep Well Block and Automated Liquid Handling Systems

The polypropylene Deep Well block is designed to be compatible with the SBS standard 96 microwell format and is therefore compatible with numerous automated liquid handling systems that accept 8 x 12 96 well assay blocks. Follow the manufacturer's recommendation for handling deep well microplates.

Examine the Drop

Carefully examine the drops under a stereo microscope (10 to 100x magnification) immediately after setting up the screen. Record all observations and be particularly careful to scan the focal plane for small crystals. Observe the drops once each day for the first week, then once a week thereafter. Records should indicate whether the drop is clear, contains precipitate, and or crystals. It is helpful to describe the drop contents using descriptive terms. Adding magnitude is also helpful. Example: 4+ yellow/brown fine precipitate, 2+ small bipyramid crystals, clear drop, 3+ needle shaped crystals in 1+ white precipitate. One may also employ a standard numerical scoring scheme (Clear = 0, Precipitate = 1, Crystal = 10, etc). Figure 6, on the left side of page 2 shows typical examples of what one might observe in a crystallization experiment.

Interpreting MembFac HT

Clear drops indicate that either the relative supersaturation of the sample and reagent is too low or the drop has not yet completed equilibration. If the drop remains clear after 3 to 4 weeks consider repeating the screen condition and doubling the sample concentration. If more than 70 of the 96 screen drops are clear consider doubling the sample concentration and repeating the entire screen.

Drops containing precipitate indicate either the relative super saturation of the sample and reagent is too high, the sample has denatured, or the sample

is heterogeneous. To reduce the relative supersaturation, dilute the sample twofold and repeat the screen condition. If more than 70 of the 96 screen drops contain precipitate and no crystals are present, consider diluting the sample concentration in half and repeating the entire screen. If sample denaturation is suspect, take measures to stabilize the sample (add reducing agent, ligands, glycerol, salt, or other stabilizing agents). If the sample is impure, aggregated, or heterogeneous take measures to pursue homogeneity. It is possible to obtain crystals from precipitate so do not discard nor ignore a drop containing precipitate. If possible, examine drops containing precipitate under polarizing optics to differentiate precipitate from microcrystalline material.

If the drop contains a macromolecular crystal the relative supersaturation of the sample and reagent is appropriate for crystal nucleation and growth. The next step is to optimize the preliminary conditions (pH, salt type, salt concentration, precipitant type, precipitant concentration, sample concentration, temperature, additives, and other crystallization variables) which produced the crystal in order to improve crystal size and quality.

Compare the observations between the 4°C and room temperature incubation to determine the effect of temperature on sample solubility. Different results in the same drops at different temperatures indicate that sample solubility is temperature dependent and that one should include temperature as a variable in subsequent screens and optimization experiments.

Retain and observe plates until the drops are dried out. Crystal growth can occur within 15 minutes or one year.

MembFac HT Formulation

Crystallization reagents are formulated using the highest purity chemicals, ultrapure water (18.2 Megohm-cm, 5 ppb TOC) and are sterile filtered using 0.22 micron filters into sterile Deep Well blocks (no preservatives added).

Crystallization reagents are readily reproduced using Hampton Research Optimize™ and StockOptions™ stock solutions of salts, polymers and buffers. Optimize and StockOptions stock reagents make reproducing crystallization screen reagents accurate, precise, fast, convenient and easy. Dilutions can be performed directly into the crystallization plate using Optimize and StockOptions stock reagents.

Crystallization reagents containing buffers are formulated by creating a 1.0 M stock buffer, titrated to the desired pH using Hydrochloric acid or Sodium hydroxide. The buffer is then diluted with the other reagent components and water. No further pH adjustment is required.

Crystallization reagents are stable at room temperature and are best if used within 12 months of receipt. To enhance reagent stability the crystallization reagents can be stored at 4°C or -20°C. Avoid ultraviolet light to preserve reagent stability.

If the sample contains phosphate, borate, or carbonate buffers it is possible to obtain inorganic crystals (false positives) when using crystallization reagents containing divalent cations such as magnesium, calcium, or zinc. To avoid false positives use phosphate, borate, or carbonate buffers at concentrations of 10 mM or less or exchange the phosphate, borate, or carbonate buffer with a more soluble buffer that does not complex with divalent cations.

References and Readings

1. Crystallization of nucleic acids and proteins, Edited by A. Ducruix and R. Giegé, The Practical Approach Series, Oxford Univ. Press, 1992.
2. Current approaches to macromolecular crystallization. McPherson, A. Eur. J. Biochem. 189, 1-23, 1990.
3. Sparse Matrix Sampling: a screening method for crystallization of proteins. Jancarik, J. and Kim, S.H. J. Appl. Cryst., 24,409-411, 1991.
4. Protein and Nucleic Acid Crystallization. Methods, A Companion to Methods in Enzymology, Academic Press, Volume 1, Number 1, August 1990.
5. Jarmila Jancarik, University of California Berkeley personal communication.
6. A comparison of salts for the crystallization of macromolecules. McPherson, A. Protein Science, 10:418-422, 2001.
7. Entering a new phase: Using solvent halide ions in protein structure determination. Dauter, Z. and Dauter, M. Structure, Vol 9, R21-26, Feb 2001.
8. Efficiency Analysis of Screening Protocols Used in Protein Crystallization, B. W. Segelke, Journal of Crystal Growth 232 : 553-562 (2001).
9. A novel approach to crystallizing proteins under oil. D'Arcy, A. et al. Journal of Crystal Growth, (1996) 168, 175-180.
10. Chayen, N. et al, J. Appl. Cryst. (1990) 23, 297.
11. Gilliland, G.L., Tung, M., Blakeslee, D.M. and Ladner, J. 1994. The Biological Macromolecule Crystallization Database, Version 3.0: New Features, Data, and the NASA Archive for Protein Crystal Growth Data. Acta Crystallogr. D50 408-413.

Technical Support

Inquiries regarding MembFac HT reagent formulation, interpretation of screen results, optimization strategies and general inquiries regarding crystallization are welcome. Please e-mail, fax, or telephone your request to Hampton Research. Fax and e-mail Technical Support are available 24 hours a day. Telephone technical support is available 8:00 a.m. to 4:30 p.m. USA Pacific Standard Time.

Hampton Research
34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

Technical Support e-mail: tech@hrmail.com

Website: www.hamptonresearch.com

| Well # | Salt | Well # | Buffer ◇ | Well # | Precipitant |
|-----------|--|-----------|--|-----------|--|
| 1. (A1) | 0.1 M Sodium chloride | 1. (A1) | 0.1 M Sodium acetate trihydrate pH 4.6 | 1. (A1) | 12% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 2. (A2) | 0.1 M Zinc acetate dihydrate | 2. (A2) | 0.1 M Sodium acetate trihydrate pH 4.6 | 2. (A2) | 12% w/v Polyethylene glycol 4,000 |
| 3. (A3) | 0.2 M Ammonium sulfate | 3. (A3) | 0.1 M Sodium acetate trihydrate pH 4.6 | 3. (A3) | 10% w/v Polyethylene glycol 4,000 |
| 4. (A4) | 0.1 M Sodium chloride | 4. (A4) | 0.1 M Sodium acetate trihydrate pH 4.6 | 4. (A4) | 12% v/v 2-Propanol |
| 5. (A5) | None | 5. (A5) | 0.1 M Sodium acetate trihydrate pH 4.6 | 5. (A5) | 12% w/v Polyethylene glycol 4,000 |
| 6. (A6) | None | 6. (A6) | 0.1 M Sodium acetate trihydrate pH 4.6 | 6. (A6) | 1.0 M Ammonium sulfate |
| 7. (A7) | None | 7. (A7) | 0.1 M Sodium acetate trihydrate pH 4.6 | 7. (A7) | 1.0 M Magnesium sulfate heptahydrate |
| 8. (A8) | 0.1 M Magnesium chloride hexahydrate | 8. (A8) | 0.1 M Sodium acetate trihydrate pH 4.6 | 8. (A8) | 18% v/v Polyethylene glycol 400 |
| 9. (A9) | 0.1 M Lithium sulfate monohydrate | 9. (A9) | 0.1 M Sodium acetate trihydrate pH 4.6 | 9. (A9) | 1.0 M Ammonium phosphate monobasic |
| 10. (A10) | 0.1 M Sodium chloride | 10. (A10) | 0.1 M Sodium acetate trihydrate pH 4.6 | 10. (A10) | 12% w/v Polyethylene glycol 6,000 |
| 11. (A11) | 0.1 M Magnesium chloride hexahydrate | 11. (A11) | 0.1 M Sodium acetate trihydrate pH 4.6 | 11. (A11) | 12% w/v Polyethylene glycol 6,000 |
| 12. (A12) | 0.1 M Sodium chloride | 12. (A12) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 12. (A12) | 18% v/v Polyethylene glycol 400 |
| 13. (B1) | 0.1 M Lithium sulfate monohydrate | 13. (B1) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 13. (B1) | 12% w/v Polyethylene glycol 4,000 |
| 14. (B2) | 0.1 M Sodium citrate tribasic dihydrate | 14. (B2) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 14. (B2) | 10% v/v 2-Propanol |
| 15. (B3) | 0.1 M Sodium chloride | 15. (B3) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 15. (B3) | 12% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 16. (B4) | None | 16. (B4) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 16. (B4) | 1.0 M Magnesium sulfate heptahydrate |
| 17. (B5) | 0.1 M Sodium chloride | 17. (B5) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 17. (B5) | 12% w/v Polyethylene glycol 4,000 |
| 18. (B6) | 0.1 M Lithium sulfate monohydrate | 18. (B6) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 18. (B6) | 12% w/v Polyethylene glycol 6,000 |
| 19. (B7) | 0.1 M Magnesium chloride hexahydrate | 19. (B7) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 19. (B7) | 4% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 20. (B8) | None | 20. (B8) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 20. (B8) | 0.1 M Sodium chloride |
| 21. (B9) | 0.1 M Lithium sulfate monohydrate | 21. (B9) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 21. (B9) | 4% v/v Polyethylene glycol 400 |
| 22. (B10) | None | 22. (B10) | 0.1 M ADA pH 6.5 | 22. (B10) | 1.0 M Ammonium sulfate |
| 23. (B11) | 0.1 M Lithium sulfate monohydrate | 23. (B11) | 0.1 M ADA pH 6.5 | 23. (B11) | 12% w/v Polyethylene glycol 4,000, 2% v/v 2-Propanol |
| 24. (B12) | None | 24. (B12) | 0.1 M ADA pH 6.5 | 24. (B12) | 1.0 M Ammonium phosphate dibasic |
| 25. (C1) | 0.1 M Magnesium chloride hexahydrate | 25. (C1) | 0.1 M ADA pH 6.5 | 25. (C1) | 12% w/v Polyethylene glycol 6,000 |
| 26. (C2) | None | 26. (C2) | 0.1 M ADA pH 6.5 | 26. (C2) | 12% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 27. (C3) | 0.1 M Lithium sulfate monohydrate | 27. (C3) | 0.1 M ADA pH 6.5 | 27. (C3) | 1.0 M Magnesium sulfate hydrate |
| 28. (C4) | 0.3 M Lithium sulfate monohydrate | 28. (C4) | 0.1 M ADA pH 6.5 | 28. (C4) | 4% v/v Polyethylene glycol 400 |
| 29. (C5) | 0.1 M Ammonium sulfate | 29. (C5) | 0.1 M HEPES sodium pH 7.5 | 29. (C5) | 0.5 M Sodium phosphate dibasic dihydrate, 0.5 M Potassium phosphate dibasic |
| 30. (C6) | 0.1 M Sodium chloride | 30. (C6) | 0.1 M HEPES sodium pH 7.5 | 30. (C6) | 10% w/v Polyethylene glycol 4,000 |
| 31. (C7) | 0.1 M Magnesium chloride hexahydrate | 31. (C7) | 0.1 M HEPES sodium pH 7.5 | 31. (C7) | 18% v/v Polyethylene glycol 400 |
| 32. (C8) | None | 32. (C8) | 0.1 M HEPES sodium pH 7.5 | 32. (C8) | 1.0 M Potassium sodium tartrate tetrahydrate |
| 33. (C9) | 0.1 M Ammonium sulfate | 33. (C9) | 0.1 M HEPES sodium pH 7.5 | 33. (C9) | 18% v/v Polyethylene glycol 400 |
| 34. (C10) | 0.1 M Ammonium sulfate | 34. (C10) | 0.1 M HEPES sodium pH 7.5 | 34. (C10) | 10% w/v Polyethylene glycol 4,000 |
| 35. (C11) | 0.1 M Sodium citrate tribasic dihydrate | 35. (C11) | 0.1 M HEPES sodium pH 7.5 | 35. (C11) | 12% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 36. (C12) | None | 36. (C12) | 0.1 M HEPES sodium pH 7.5 | 36. (C12) | 1.0 M Sodium citrate tribasic dihydrate |
| 37. (D1) | 0.6 M Magnesium sulfate hydrate | 37. (D1) | 0.1 M HEPES sodium pH 7.5 | 37. (D1) | 4% v/v Polyethylene glycol 400 |
| 38. (D2) | 0.6 M Magnesium sulfate hydrate | 38. (D2) | 0.1 M HEPES sodium pH 7.5 | 38. (D2) | 4% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 39. (D3) | 0.1 M Lithium sulfate monohydrate | 39. (D3) | 0.1 M HEPES sodium pH 7.5 | 39. (D3) | 0.1 M Potassium sodium tartrate tetrahydrate |
| 40. (D4) | 0.1 M Lithium sulfate monohydrate | 40. (D4) | 0.1 M TRIS hydrochloride pH 8.5 | 40. (D4) | 12% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 41. (D5) | 0.1 M Ammonium phosphate dibasic | 41. (D5) | 0.1 M TRIS hydrochloride pH 8.5 | 41. (D5) | 0.5 M Sodium phosphate dibasic dihydrate, 0.5 M Potassium phosphate dibasic |
| 42. (D6) | None | 42. (D6) | 0.1 M TRIS hydrochloride pH 8.5 | 42. (D6) | 0.1 M Sodium acetate trihydrate |
| 43. (D7) | None | 43. (D7) | 0.1 M TRIS hydrochloride pH 8.5 | 43. (D7) | 0.1 M Sodium chloride |
| 44. (D8) | 0.1 M Ammonium phosphate dibasic | 44. (D8) | 0.1 M TRIS hydrochloride pH 8.5 | 44. (D8) | 12% w/v Polyethylene glycol 6,000 |
| 45. (D9) | 0.1 M Potassium sodium tartrate tetrahydrate | 45. (D9) | 0.1 M TRIS hydrochloride pH 8.5 | 45. (D9) | 0.4 M Magnesium sulfate hydrate |
| 46. (D10) | None | 46. (D10) | 0.1 M TRIS hydrochloride pH 8.5 | 46. (D10) | 0.2 M Lithium sulfate monohydrate |
| 47. (D11) | None | 47. (D11) | 0.1 M TRIS hydrochloride pH 8.5 | 47. (D11) | 0.5 M Ammonium sulfate |
| 48. (D12) | 0.1 M Sodium citrate tribasic dihydrate | 48. (D12) | 0.1 M TRIS hydrochloride pH 8.5 | 48. (D12) | 5% v/v Polyethylene glycol 400 |

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

*MembFac HT™ (Deep Well Block) contains ninety-six unique reagents beginning at position A1.
To determine the formulation of each reagent, simply read across the page.*

| Well # | Salt | Well # | Buffer ◇ | Well # | Precipitant |
|-----------|---|-----------|--|-----------|--|
| 49. (E1) | 0.02 M Calcium chloride dihydrate | 49. (E1) | 0.1 M Sodium acetate trihydrate pH 4.6 | 49. (E1) | 15% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 50. (E2) | None | 50. (E2) | None | 50. (E2) | 0.2 M Potassium sodium tartrate tetrahydrate |
| 51. (E3) | None | 51. (E3) | None | 51. (E3) | 0.2 M Ammonium phosphate monobasic |
| 52. (E4) | None | 52. (E4) | 0.1 M TRIS hydrochloride pH 8.5 | 52. (E4) | 1.0 M Ammonium sulfate |
| 53. (E5) | 0.2 M Sodium citrate tribasic dihydrate | 53. (E5) | 0.1 M HEPES sodium pH 7.5 | 53. (E5) | 15% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 54. (E6) | 0.2 M Magnesium chloride hexahydrate | 54. (E6) | 0.1 M TRIS hydrochloride pH 8.5 | 54. (E6) | 15% w/v Polyethylene glycol 4,000 |
| 55. (E7) | None | 55. (E7) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 55. (E7) | 0.7 M Sodium acetate trihydrate |
| 56. (E8) | 0.2 M Sodium citrate tribasic dihydrate | 56. (E8) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 56. (E8) | 15% v/v 2-Propanol |
| 57. (E9) | 0.2 M Ammonium acetate | 57. (E9) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 57. (E9) | 15% w/v Polyethylene glycol 4,000 |
| 58. (E10) | 0.2 M Ammonium acetate | 58. (E10) | 0.1 M Sodium acetate trihydrate pH 4.6 | 58. (E10) | 15% w/v Polyethylene glycol 4,000 |
| 59. (E11) | None | 59. (E11) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 59. (E11) | 0.5 M Ammonium phosphate monobasic |
| 60. (E12) | 0.2 M Magnesium chloride hexahydrate | 60. (E12) | 0.1 M HEPES sodium pH 7.5 | 60. (E12) | 15% v/v 2-Propanol |
| 61. (F1) | 0.2 M Sodium citrate tribasic dihydrate | 61. (F1) | 0.1 M TRIS hydrochloride pH 8.5 | 61. (F1) | 15% v/v Polyethylene glycol 400 |
| 62. (F2) | 0.2 M Calcium chloride dihydrate | 62. (F2) | 0.1 M HEPES sodium pH 7.5 | 62. (F2) | 14% v/v Polyethylene glycol 400 |
| 63. (F3) | 0.2 M Ammonium sulfate | 63. (F3) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 63. (F3) | 15% w/v Polyethylene glycol 8,000 |
| 64. (F4) | None | 64. (F4) | 0.1 M HEPES sodium pH 7.5 | 64. (F4) | 0.75 M Lithium sulfate monohydrate |
| 65. (F5) | 0.2 M Lithium sulfate monohydrate | 65. (F5) | 0.1 M TRIS hydrochloride pH 8.5 | 65. (F5) | 15% w/v Polyethylene glycol 4,000 |
| 66. (F6) | 0.2 M Magnesium acetate tetrahydrate | 66. (F6) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 66. (F6) | 10% w/v Polyethylene glycol 8,000 |
| 67. (F7) | 0.2 M Ammonium acetate | 67. (F7) | 0.1 M TRIS hydrochloride pH 8.5 | 67. (F7) | 15% v/v 2-Propanol |
| 68. (F8) | 0.2 M Ammonium sulfate | 68. (F8) | 0.1 M Sodium acetate trihydrate pH 4.6 | 68. (F8) | 12.5% w/v Polyethylene glycol 4,000 |
| 69. (F9) | 0.2 M Magnesium acetate tetrahydrate | 69. (F9) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 69. (F9) | 15% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 70. (F10) | 0.2 M Sodium acetate trihydrate | 70. (F10) | 0.1 M TRIS hydrochloride pH 8.5 | 70. (F10) | 15% w/v Polyethylene glycol 4,000 |
| 71. (F11) | 0.2 M Magnesium chloride hexahydrate | 71. (F11) | 0.1 M HEPES sodium pH 7.5 | 71. (F11) | 15% v/v Polyethylene glycol 400 |
| 72. (F12) | 0.2 M Calcium chloride dihydrate | 72. (F12) | 0.1 M Sodium acetate trihydrate pH 4.6 | 72. (F12) | 10% v/v 2-Propanol |
| 73. (G1) | None | 73. (G1) | 0.1 M Imidazole pH 6.5 | 73. (G1) | 0.5 M Sodium acetate trihydrate |
| 74. (G2) | 0.2 M Ammonium acetate | 74. (G2) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 74. (G2) | 15% v/v (+/-)-2-Methyl-2,4-pentanediol |
| 75. (G3) | 0.2 M Sodium citrate tribasic dihydrate | 75. (G3) | 0.1 M HEPES sodium pH 7.5 | 75. (G3) | 10% v/v 2-Propanol |
| 76. (G4) | 0.2 M Sodium acetate trihydrate | 76. (G4) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 76. (G4) | 15% w/v Polyethylene glycol 8,000 |
| 77. (G5) | None | 77. (G5) | 0.1 M HEPES sodium pH 7.5 | 77. (G5) | 0.4 M Potassium sodium tartrate tetrahydrate |
| 78. (G6) | 0.2 M Ammonium sulfate | 78. (G6) | None | 78. (G6) | 15% w/v Polyethylene glycol 8,000 |
| 79. (G7) | 0.2 M Ammonium sulfate | 79. (G7) | None | 79. (G7) | 15% w/v Polyethylene glycol 4,000 |
| 80. (G8) | None | 80. (G8) | None | 80. (G8) | 1.0 M Ammonium sulfate |
| 81. (G9) | None | 81. (G9) | None | 81. (G9) | 2.0 M Sodium formate |
| 82. (G10) | None | 82. (G10) | 0.1 M Sodium acetate trihydrate pH 4.6 | 82. (G10) | 1.0 M Sodium formate |
| 83. (G11) | None | 83. (G11) | 0.1 M HEPES sodium pH 7.5 | 83. (G11) | 0.4 M Sodium phosphate monobasic monohydrate, 0.4 M Potassium phosphate monobasic |
| 84. (G12) | None | 84. (G12) | 0.1 M TRIS hydrochloride pH 8.5 | 84. (G12) | 4% w/v Polyethylene glycol 8,000 |
| 85. (H1) | None | 85. (H1) | 0.1 M Sodium acetate trihydrate pH 4.6 | 85. (H1) | 4% w/v Polyethylene glycol 4,000 |
| 86. (H2) | None | 86. (H2) | 0.1 M HEPES sodium pH 7.5 | 86. (H2) | 0.7 M Sodium citrate tribasic dihydrate |
| 87. (H3) | None | 87. (H3) | 0.1 M HEPES sodium pH 7.5 | 87. (H3) | 2% v/v Polyethylene glycol 400, 1.0 M Ammonium sulfate |
| 88. (H4) | None | 88. (H4) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6 | 88. (H4) | 10% v/v 2-Propanol, 10% w/v Polyethylene glycol 4,000 |
| 89. (H5) | None | 89. (H5) | 0.1 M HEPES sodium pH 7.5 | 89. (H5) | 5% v/v 2-Propanol, 10% w/v Polyethylene glycol 4,000 |
| 90. (H6) | 0.05 M Potassium phosphate monobasic | 90. (H6) | None | 90. (H6) | 10% w/v Polyethylene glycol 8,000 |
| 91. (H7) | None | 91. (H7) | None | 91. (H7) | 15% w/v Polyethylene glycol 1,500 |
| 92. (H8) | None | 92. (H8) | None | 92. (H8) | 0.1 M Magnesium formate dihydrate |
| 93. (H9) | 0.2 M Zinc acetate dihydrate | 93. (H9) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 93. (H9) | 9% w/v Polyethylene glycol 8,000 |
| 94. (H10) | 0.2 M Calcium acetate hydrate | 94. (H10) | 0.1 M Sodium cacodylate trihydrate pH 6.5 | 94. (H10) | 9% w/v Polyethylene glycol 8,000 |
| 95. (H11) | None | 95. (H11) | 0.1 M Sodium acetate trihydrate pH 4.6 | 95. (H11) | 1.0 M Ammonium sulfate |
| 96. (H12) | None | 96. (H12) | 0.1 M TRIS hydrochloride pH 8.5 | 96. (H12) | 1.0 M Ammonium phosphate monobasic |

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

*MembFac HT™ (Deep Well Block) contains ninety-six unique reagents beginning at position A1.
To determine the formulation of each reagent, simply read across the page.*

**HAMPTON
RESEARCH**

Solutions for Crystal Growth

Sample: _____ Sample Concentration: _____
 Sample Buffer: _____ Date: _____
 Reservoir Volume: _____ Temperature: _____
 Drop Volume: Total _____ μ l Sample _____ μ l Reservoir _____ μ l Additive _____ μ l

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals
- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2 mm)
- 9 Single Crystals (3D Growth > 0.2 mm)

MembFac HT™ - HR2-137 Scoring Sheet

Date: Date: Date:

| | | | | |
|-----------|--|--|--|--|
| 1. (A1) | 0.1 M Sodium chloride, 0.1 M Sodium acetate trihydrate pH 4.6, 12% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 2. (A2) | 0.1 M Zinc acetate dihydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 12% w/v Polyethylene glycol 4,000 | | | |
| 3. (A3) | 0.2 M Ammonium sulfate, 0.1 M Sodium acetate trihydrate pH 4.6, 10% w/v Polyethylene glycol 4,000 | | | |
| 4. (A4) | 0.1 M Sodium chloride, 0.1 M Sodium acetate trihydrate pH 4.6, 12% v/v 2-Propanol | | | |
| 5. (A5) | 0.1 M Sodium acetate trihydrate pH 4.6, 12% w/v Polyethylene glycol 4,000 | | | |
| 6. (A6) | 0.1 M Sodium acetate trihydrate pH 4.6, 1.0 M Ammonium sulfate | | | |
| 7. (A7) | 0.1 M Sodium acetate trihydrate pH 4.6, 1.0 M Magnesium sulfate heptahydrate | | | |
| 8. (A8) | 0.1 M Magnesium chloride hexahydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 18% v/v Polyethylene glycol 400 | | | |
| 9. (A9) | 0.1 M Lithium sulfate monohydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 1.0 M Ammonium phosphate monobasic | | | |
| 10. (A10) | 0.1 M Sodium chloride, 0.1 M Sodium acetate trihydrate pH 4.6, 12% w/v Polyethylene glycol 6,000 | | | |
| 11. (A11) | 0.1 M Magnesium chloride hexahydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 12% w/v Polyethylene glycol 6,000 | | | |
| 12. (A12) | 0.1 M Sodium chloride, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 18% v/v Polyethylene glycol 400 | | | |
| 13. (B1) | 0.1 M Lithium sulfate monohydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 12% w/v Polyethylene glycol 4,000 | | | |
| 14. (B2) | 0.1 M Sodium citrate tribasic dihydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 10% v/v 2-Propanol | | | |
| 15. (B3) | 0.1 M Sodium chloride, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 12% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 16. (B4) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 1.0 M Magnesium sulfate heptahydrate | | | |
| 17. (B5) | 0.1 M Sodium chloride, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 12% w/v Polyethylene glycol 4,000 | | | |
| 18. (B6) | 0.1 M Lithium sulfate monohydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 12% w/v Polyethylene glycol 6,000 | | | |
| 19. (B7) | 0.1 M Magnesium chloride hexahydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 4% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 20. (B8) | 0.1 M Sodium citrate trihydrate dihydrate pH 5.6, 0.1 M Sodium chloride | | | |
| 21. (B9) | 0.1 M Lithium sulfate monohydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 4% v/v Polyethylene glycol 400 | | | |
| 22. (B10) | 0.1 M ADA pH 6.5, 1.0 M Ammonium sulfate | | | |
| 23. (B11) | 0.1 M Lithium sulfate monohydrate, 0.1 M ADA pH 6.5, 12% w/v Polyethylene glycol 4,000, 2% v/v 2-Propanol | | | |
| 24. (B12) | 0.1 M ADA pH 6.5, 1.0 M Ammonium phosphate dibasic | | | |
| 25. (C1) | 0.1 M Magnesium chloride hexahydrate, 0.1 M ADA pH 6.5, 12% w/v Polyethylene glycol 6,000 | | | |
| 26. (C2) | 0.1 M ADA pH 6.5, 12% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 27. (C3) | 0.1 M Lithium sulfate monohydrate, 0.1 M ADA pH 6.5, 1.0 M Magnesium sulfate hydrate | | | |
| 28. (C4) | 0.3 M Lithium sulfate monohydrate, 0.1 M ADA pH 6.5, 4% v/v Polyethylene glycol 400 | | | |
| 29. (C5) | 0.1 M Ammonium sulfate, 0.1 M HEPES sodium pH 7.5, 0.5 M Sodium phosphate dibasic dihydrate, 0.5 M Potassium phosphate dibasic | | | |
| 30. (C6) | 0.1 M Sodium chloride, 0.1 M HEPES sodium pH 7.5, 10% w/v Polyethylene glycol 4,000 | | | |
| 31. (C7) | 0.1 M Magnesium chloride hexahydrate, 0.1 M HEPES sodium pH 7.5, 18% v/v Polyethylene glycol 400 | | | |
| 32. (C8) | 0.1 M HEPES sodium pH 7.5, 1.0 M Potassium sodium tartrate tetrahydrate | | | |
| 33. (C9) | 0.1 M Ammonium sulfate, 0.1 M HEPES sodium pH 7.5, 18% v/v Polyethylene glycol 400 | | | |
| 34. (C10) | 0.1 M Ammonium sulfate, 0.1 M HEPES sodium pH 7.5, 10% w/v Polyethylene glycol 4,000 | | | |
| 35. (C11) | 0.1 M Sodium citrate tribasic dihydrate, 0.1 M HEPES sodium pH 7.5, 12% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 36. (C12) | 0.1 M HEPES sodium pH 7.5, 1.0 M Sodium citrate tribasic dihydrate | | | |
| 37. (D1) | 0.6 M Magnesium sulfate hydrate, 0.1 M HEPES sodium pH 7.5, 4% v/v Polyethylene glycol 400 | | | |
| 38. (D2) | 0.6 M Magnesium sulfate hydrate, 0.1 M HEPES sodium pH 7.5, 4% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 39. (D3) | 0.1 M Lithium sulfate monohydrate, 0.1 M HEPES sodium pH 7.5, 0.1 M Potassium sodium tartrate tetrahydrate | | | |
| 40. (D4) | 0.1 M Lithium sulfate monohydrate, 0.1 M TRIS hydrochloride pH 8.5, 12% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 41. (D5) | 0.1 M Ammonium phosphate dibasic, 0.1 M TRIS hydrochloride pH 8.5, 0.5 M Na phosphate dibasic dihydrate, 0.5 M Potassium phosphate dibasic | | | |
| 42. (D6) | 0.1 M TRIS hydrochloride pH 8.5, 0.1 M Sodium acetate trihydrate | | | |
| 43. (D7) | 0.1 M TRIS hydrochloride pH 8.5, 0.1 M Sodium chloride | | | |
| 44. (D8) | 0.1 M Ammonium phosphate dibasic, 0.1 M TRIS hydrochloride pH 8.5, 12% w/v Polyethylene glycol 6,000 | | | |
| 45. (D9) | 0.1 M Potassium sodium tartrate tetrahydrate, 0.1 M TRIS hydrochloride pH 8.5, 0.4 M Magnesium sulfate hydrate | | | |
| 46. (D10) | 0.1 M TRIS hydrochloride pH 8.5, 0.2 M Lithium sulfate monohydrate | | | |
| 47. (D11) | 0.1 M TRIS hydrochloride pH 8.5, 0.5 M Ammonium sulfate | | | |
| 48. (D12) | 0.1 M Sodium citrate tribasic dihydrate, 0.1 M TRIS hydrochloride pH 8.5, 5% v/v Polyethylene glycol 400 | | | |

Sample: _____ Sample Concentration: _____
 Sample Buffer: _____ Date: _____
 Reservoir Volume: _____ Temperature: _____
 Drop Volume: Total _____ μ l Sample _____ μ l Reservoir _____ μ l Additive _____ μ l

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals
- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2mm)
- 9 Single Crystals (3D Growth > 0.2mm)

MembFac HT™ - HR2-137 Scoring Sheet

Date: Date: Date:

| | | | | |
|-----------|---|--|--|--|
| 49. (E1) | 0.02 M Calcium chloride dihydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 15% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 50. (E2) | 0.2 M Potassium sodium tartrate tetrahydrate | | | |
| 51. (E3) | 0.2 M Ammonium phosphate monobasic | | | |
| 52. (E4) | 0.1 M TRIS hydrochloride pH 8.5, 1.0 M Ammonium sulfate | | | |
| 53. (E5) | 0.2 M Sodium citrate tribasic dihydrate, 0.1 M HEPES sodium pH 7.5, 15% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 54. (E6) | 0.2 M Magnesium chloride hexahydrate, 0.1 M TRIS hydrochloride pH 8.5, 15% w/v Polyethylene glycol 4,000 | | | |
| 55. (E7) | 0.1 M Sodium cacodylate trihydrate pH 6.5, 0.7 M Sodium acetate trihydrate | | | |
| 56. (E8) | 0.2 M Sodium citrate tribasic dihydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 15% v/v 2-Propanol | | | |
| 57. (E9) | 0.2 M Ammonium acetate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 15% w/v Polyethylene glycol 4,000 | | | |
| 58. (E10) | 0.2 M Ammonium acetate, 0.1 M Sodium acetate trihydrate pH 4.6, 15% w/v Polyethylene glycol 4,000 | | | |
| 59. (E11) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 0.5 M Ammonium phosphate monobasic | | | |
| 60. (E12) | 0.2 M Magnesium chloride hexahydrate, 0.1 M HEPES sodium pH 7.5, 15% v/v 2-Propanol | | | |
| 61. (F1) | 0.2 M Sodium citrate tribasic dihydrate, 0.1 M TRIS hydrochloride pH 8.5, 15% v/v Polyethylene glycol 400 | | | |
| 62. (F2) | 0.2 M Calcium chloride dihydrate, 0.1 M HEPES sodium pH 7.5, 14% v/v Polyethylene glycol 400 | | | |
| 63. (F3) | 0.2 M Ammonium sulfate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 15% w/v Polyethylene glycol 8,000 | | | |
| 64. (F4) | 0.1 M HEPES sodium pH 7.5, 0.75 M Lithium sulfate monohydrate | | | |
| 65. (F5) | 0.2 M Lithium sulfate monohydrate, 0.1 M TRIS hydrochloride pH 8.5, 15% w/v Polyethylene glycol 4,000 | | | |
| 66. (F6) | 0.2 M Magnesium acetate tetrahydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 10% w/v Polyethylene glycol 8,000 | | | |
| 67. (F7) | 0.2 M Ammonium acetate, 0.1 M TRIS hydrochloride pH 8.5, 15% v/v 2-Propanol | | | |
| 68. (F8) | 0.2 M Ammonium sulfate, 0.1 M Sodium acetate trihydrate pH 4.6, 12.5% w/v Polyethylene glycol 4,000 | | | |
| 69. (F9) | 0.2 M Magnesium acetate tetrahydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 15% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 70. (F10) | 0.2 M Sodium acetate trihydrate, 0.1 M TRIS hydrochloride pH 8.5, 15% w/v Polyethylene glycol 4,000 | | | |
| 71. (F11) | 0.2 M Magnesium chloride hexahydrate, 0.1 M HEPES sodium pH 7.5, 15% v/v Polyethylene glycol 400 | | | |
| 72. (F12) | 0.2 M Calcium chloride dihydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 10% v/v 2-Propanol | | | |
| 73. (G1) | 0.1 M Imidazole pH 6.5, 0.5 M Sodium acetate trihydrate | | | |
| 74. (G2) | 0.2 M Ammonium acetate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 15% v/v (+/-)-2-Methyl-2,4-pentanediol | | | |
| 75. (G3) | 0.2 M Sodium citrate tribasic dihydrate, 0.1 M HEPES sodium pH 7.5, 10% v/v 2-Propanol | | | |
| 76. (G4) | 0.2 M Sodium acetate trihydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 15% w/v Polyethylene glycol 8,000 | | | |
| 77. (G5) | 0.1 M HEPES sodium pH 7.5, 0.4 M Potassium sodium tartrate tetrahydrate | | | |
| 78. (G6) | 0.2 M Ammonium sulfate, 15% w/v Polyethylene glycol 8,000 | | | |
| 79. (G7) | 0.2 M Ammonium sulfate, 15% w/v Polyethylene glycol 4,000 | | | |
| 80. (G8) | 1.0 M Ammonium sulfate | | | |
| 81. (G9) | 2.0 M Sodium formate | | | |
| 82. (G10) | 0.1 M Sodium acetate trihydrate pH 4.6, 1.0 M Sodium formate | | | |
| 83. (G11) | 0.1 M HEPES sodium pH 7.5, 0.4 M Sodium phosphate monobasic monohydrate, 0.4 M Potassium phosphate monobasic | | | |
| 84. (G12) | 0.1 M TRIS hydrochloride pH 8.5, 4% w/v Polyethylene glycol 8,000 | | | |
| 85. (H1) | 0.1 M Sodium acetate trihydrate pH 4.6, 4% w/v Polyethylene glycol 4,000 | | | |
| 86. (H2) | 0.1 M HEPES sodium pH 7.5, 0.7 M Sodium citrate tribasic dihydrate | | | |
| 87. (H3) | 0.1 M HEPES sodium pH 7.5, 2% v/v Polyethylene glycol 400, 1.0 M Ammonium sulfate | | | |
| 88. (H4) | 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 10% v/v 2-Propanol, 10% w/v Polyethylene glycol 4,000 | | | |
| 89. (H5) | 0.1 M HEPES sodium pH 7.5, 5% v/v 2-Propanol, 10% w/v Polyethylene glycol 4,000 | | | |
| 90. (H6) | 0.05 M Potassium phosphate monobasic, 10% w/v Polyethylene glycol 8,000 | | | |
| 91. (H7) | 15% w/v Polyethylene glycol 1,500 | | | |
| 92. (H8) | 0.1 M Magnesium formate dihydrate | | | |
| 93. (H9) | 0.2 M Zinc acetate dihydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 9% w/v Polyethylene glycol 8,000 | | | |
| 94. (H10) | 0.2 M Calcium acetate hydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 9% w/v Polyethylene glycol 8,000 | | | |
| 95. (H11) | 0.1 M Sodium acetate trihydrate pH 4.6, 1.0 M Ammonium sulfate | | | |
| 96. (H12) | 0.1 M TRIS hydrochloride pH 8.5, 1.0 M Ammonium phosphate monobasic | | | |