


Additive Screen™

Additive Screen™ is a kit designed to allow the rapid and convenient evaluation of 96 unique additives and their ability to influence the crystallization of the sample. The screen is designed to be compatible with most popular crystallization reagents including all reagents utilized in all of the Hampton Research screens.

Each of the additives is preformulated in deionized water and sterile filtered using a 0.2 micron filter. Recommended storage for the Additive Screen kit is -20 to 4°C. Allow the kit to equilibrate to room temperature prior to removing the cap from the tube. If reagents precipitate during cold storage, warm the tube at 37°C for up to 60 minutes and invert several times to solubilize the reagents.

1




The Additive Screen™ kit is a complete reagent kit designed to provide a rapid screening method for the manipulation of sample-sample and sample solvent interactions to enhance or alter sample solubility.

The Additive Screen evaluates the manipulation factors of multivalent cations, salts, amino acid, dissociating agents, linkers, polyamines, chaotropes, co-factors, reducing agents, polymers, chelating agent, carbohydrates, polyols, non-detergents, amphiphiles, detergents, osmolyte, organic (non-volatile) and organic (volatile) reagents.

The Additive Screen kit is to be used before and during the optimization of preliminary crystallization conditions.

Each Additive Screen kit contains 1 milliliter of 96 unique additives formulated to allow one to rapidly screen with less than 100 microliters of sample.

2



This guide will describe the use of the Additive Screen kit using the Sitting Drop Vapor Diffusion method and a 1 milliliter reservoir volume. Other methods such as Hanging Drop Vapor Diffusion crystallization, and MicroBatch may also be utilized as well as smaller reservoir and drop volumes. A complete description of the Hanging, Sitting, Sandwich Drop, Dialysis and other crystallization methods are available from the Hampton Research Crystal Growth 101 Library.

Two separate methods of setup are to be used for volatile and non-volatile additives.

Reservoir setup for **non-volatile** Additives (A1 - G8):


A. Pipet 1 milliliter of crystallization reagent into the reservoir only.

Reservoir setup for **volatile** Additives (G9 - H12):

A. Pipet 900 µl of crystallization reagent into the reservoir.

B. Pipet and mix 100 µl of the volatile additive into the reservoir.

3



Drop setup for **non-volatile** Additives

A. Pipet 5 µl of sample onto a sitting drop post.

B. Pipet 1 µl of additive into the sample drop.

C. Pipet 4 µl of the crystallization reagent into the sample/additive drop.

D. Seal the reservoir with tape or grease and slides.

E. Repeat for remaining additives.

Drop setup for **volatile** Additives

A. Pipet 5 µl of sample onto a sitting drop post.


B. Pipet 5 µl of the crystallization reagent/additive mixture from the reservoir into the sample drop.

C. Seal the reservoir with tape or grease and slides.


D. Repeat for the remaining additives.

Figure 1


Typical observations in a crystallization experiment




Clear Drop




Skin/Precipitate




Precipitate




Precipitate/Phase




Quasi Crystals



Microcrystals




Needle Cluster



Plates



Rod Cluster



Single Crystal

4

Examine The Drop

Carefully examine the drops under a stereo microscope (10 to 100x magnification) immediately after setting up the screen. Record all observations and be particularly careful to scan the focal plane for small crystals. Observe the drops once each day for the first week, then once a week thereafter. Records should indicate whether the drop is clear, contains precipitate, and or crystals. It is helpful to describe the drop contents using descriptive terms. Adding magnitude is also helpful. Example: 4+ yellow/brown fine precipitate, 2+ small bipyramid crystals, clear drop, 3+ needle shaped crystals

5

References and Readings

1. Crystallization of membrane proteins. Edited by Hartmut Michel, CRC Press, 1991.
2. Crystallization of nucleic acids and proteins, Edited by A. Ducruix and R. Giege, The Practical Approach Series, Oxford Univ. Press, 1992 175-191.
3. Screening and optimization strategies for macromolecular crystal growth. Cudney, B. et al, Acta Cryst. (1994). D50, 414-423.
4. Use of glycerol, polyols and other protein structure stabilizing agents in protein crystallization. R. Sousa. Acta Cryst. (1995) D51, 271-277.
5. Influence of divalent cations on protein crystallization. Trakhanov, S. and Quioco, F.A. (1995) Protein Science 4(9): 1914-1919.
6. Non-detergent sulphobetaines: a new class of mild solubilizing agents for protein purification. L. Vuillard, C. Braun-Breton, T. Rabilloud, Biochem. J. (1995) 305, 337-343.
7. A new additive for protein crystallization. L. Vuillard, T. Rabilloud, R. Leberman, C. Berthet-Colominas, St. Cusack. FEBS Letters, 353 (1994) 294-296.

in 1+ white precipitate. One may also employ a standard numerical scoring scheme (Clear = 0, Precipitate = 1, Crystal = 10, etc). Figure 1 (left side of page 2) shows typical examples of what one might observe in a crystallization experiment.

Technical Support

Inquiries regarding Additive Screen reagent formulation, interpretation of screen results, optimization strategies and general inquiries regarding crystallization are welcome. Please e-mail, fax, or telephone your request to Hampton Research. Fax and e-mail Technical Support are available 24 hours a day. Telephone technical support is available 8:00 a.m. to 4:30 p.m. USA Pacific Standard Time.

Hampton Research
34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

Technical Support e-mail: tech@hrmail.com

Website: www.hamptonresearch.com

| Tube # | Salt | Tube # | Classification | Tube # | Suggested Drop Concentration |
|-----------|---|-----------|----------------------|-----------|------------------------------|
| 1. (A1) | 0.1 M Barium chloride dihydrate | 1. (A1) | Multivalent | 1. (A1) | 0.01 M (10 mM) |
| 2. (A2) | 0.1 M Cadmium chloride hydrate | 2. (A2) | Multivalent | 2. (A2) | 0.01 M (10 mM) |
| 3. (A3) | 0.1 M Calcium chloride dihydrate | 3. (A3) | Multivalent | 3. (A3) | 0.01 M (10 mM) |
| 4. (A4) | 0.1 M Cobalt(II) chloride hexahydrate | 4. (A4) | Multivalent | 4. (A4) | 0.01 M (10 mM) |
| 5. (A5) | 0.1 M Copper(II) chloride dihydrate | 5. (A5) | Multivalent | 5. (A5) | 0.01 M (10 mM) |
| 6. (A6) | 0.1 M Magnesium chloride hexahydrate | 6. (A6) | Multivalent | 6. (A6) | 0.01 M (10 mM) |
| 7. (A7) | 0.1 M Manganese(II) chloride tetrahydrate | 7. (A7) | Multivalent | 7. (A7) | 0.01 M (10 mM) |
| 8. (A8) | 0.1 M Strontium chloride hexahydrate | 8. (A8) | Multivalent | 8. (A8) | 0.01 M (10 mM) |
| 9. (A9) | 0.1 M Yttrium(III) chloride hexahydrate | 9. (A9) | Multivalent | 9. (A9) | 0.01 M (10 mM) |
| 10. (A10) | 0.1 M Zinc chloride | 10. (A10) | Multivalent | 10. (A10) | 0.01 M (10 mM) |
| 11. (A11) | 0.1 M Iron(III) chloride hexahydrate | 11. (A11) | Multivalent | 11. (A11) | 0.01 M (10 mM) |
| 12. (A12) | 0.1 M Nickel(II) chloride hexahydrate | 12. (A12) | Multivalent | 12. (A12) | 0.01 M (10 mM) |
| 13. (B1) | 0.1 M Chromium(III) chloride hexahydrate | 13. (B1) | Multivalent | 13. (B1) | 0.01 M (10 mM) |
| 14. (B2) | 0.1 M Praseodymium(III) acetate hydrate | 14. (B2) | Multivalent | 14. (B2) | 0.01 M (10 mM) |
| 15. (B3) | 1.0 M Ammonium sulfate | 15. (B3) | Salt | 15. (B3) | 0.1 M (100 mM) |
| 16. (B4) | 1.0 M Potassium chloride | 16. (B4) | Salt | 16. (B4) | 0.1 M (100 mM) |
| 17. (B5) | 1.0 M Lithium chloride | 17. (B5) | Salt | 17. (B5) | 0.1 M (100 mM) |
| 18. (B6) | 2.0 M Sodium chloride | 18. (B6) | Salt | 18. (B6) | 0.2 M (200 mM) |
| 19. (B7) | 0.5 M Sodium fluoride | 19. (B7) | Salt | 19. (B7) | 0.05 M (50 mM) |
| 20. (B8) | 1.0 M Sodium iodide | 20. (B8) | Salt | 20. (B8) | 0.1 M (100 mM) |
| 21. (B9) | 2.0 M Sodium thiocyanate | 21. (B9) | Salt | 21. (B9) | 0.2 M (200 mM) |
| 22. (B10) | 1.0 M Potassium sodium tartrate tetrahydrate | 22. (B10) | Salt | 22. (B10) | 0.1 M (100 mM) |
| 23. (B11) | 1.0 M Sodium citrate tribasic dihydrate | 23. (B11) | Salt | 23. (B11) | 0.1 M (100 mM) |
| 24. (B12) | 1.0 M Cesium chloride | 24. (B12) | Salt | 24. (B12) | 0.1 M (100 mM) |
| 25. (C1) | 1.0 M Sodium malonate pH 7.0 | 25. (C1) | Salt | 25. (C1) | 0.1 M (100 mM) |
| 26. (C2) | 0.1 M L-Proline | 26. (C2) | Amino Acid | 26. (C2) | 0.01 M (10 mM) |
| 27. (C3) | 0.1 M Phenol | 27. (C3) | Dissociating Agent | 27. (C3) | 0.01 M (10 mM) |
| 28. (C4) | 30% v/v Dimethyl sulfoxide | 28. (C4) | Dissociating Agent | 28. (C4) | 3.0% |
| 29. (C5) | 0.1 M Sodium bromide | 29. (C5) | Dissociating Agent | 29. (C5) | 0.01 M (10 mM) |
| 30. (C6) | 30% w/v 6-Aminohexanoic acid | 30. (C6) | Linker | 30. (C6) | 3.0% |
| 31. (C7) | 30% w/v 1,5-Diaminopentane dihydrochloride | 31. (C7) | Linker | 31. (C7) | 3.0% |
| 32. (C8) | 30% w/v 1,6-Diaminohexane | 32. (C8) | Linker | 32. (C8) | 3.0% |
| 33. (C9) | 30% w/v 1,8-Diaminooctane | 33. (C9) | Linker | 33. (C9) | 3.0% |
| 34. (C10) | 1.0 M Glycine | 34. (C10) | Linker | 34. (C10) | 0.1 M (100 mM) |
| 35. (C11) | 0.3 M Glycyl-glycyl-glycine | 35. (C11) | Linker | 35. (C11) | 0.03 M (30 mM) |
| 36. (C12) | 0.1 M Taurine | 36. (C12) | Linker | 36. (C12) | 0.01 M (10 mM) |
| 37. (D1) | 0.1 M Betaine hydrochloride | 37. (D1) | Linker | 37. (D1) | 0.01 M (10 mM) |
| 38. (D2) | 0.1 M Spermidine | 38. (D2) | Polyamine | 38. (D2) | 0.01 M (10 mM) |
| 39. (D3) | 0.1 M Spermine tetrahydrochloride | 39. (D3) | Polyamine | 39. (D3) | 0.01 M (10 mM) |
| 40. (D4) | 0.1 M Hexamine cobalt(III) chloride | 40. (D4) | Polyamine | 40. (D4) | 0.01 M (10 mM) |
| 41. (D5) | 0.1 M Sarcosine | 41. (D5) | Polyamine / Osmolyte | 41. (D5) | 0.01 M (10 mM) |
| 42. (D6) | 0.1 M Trimethylamine hydrochloride | 42. (D6) | Chaotrope | 42. (D6) | 0.01 M (10 mM) |
| 43. (D7) | 1.0 M Guanidine hydrochloride | 43. (D7) | Chaotrope | 43. (D7) | 0.1 M (100 mM) |
| 44. (D8) | 0.1 M Urea | 44. (D8) | Chaotrope | 44. (D8) | 0.01 M (10 mM) |
| 45. (D9) | 0.1 M β -Nicotinamide adenine dinucleotide hydrate | 45. (D9) | Co-factor | 45. (D9) | 0.01 M (10 mM) |
| 46. (D10) | 0.1 M Adenosine-5'-triphosphate disodium salt hydrate | 46. (D10) | Co-factor | 46. (D10) | 0.01 M (10 mM) |
| 47. (D11) | 0.1 M TCEP hydrochloride | 47. (D11) | Reducing Agent | 47. (D11) | 0.01 M (10 mM) |
| 48. (D12) | 0.01 M GSH (L-Glutathione reduced), 0.01 M GSSG (L-Glutathione oxidized) | 48. (D12) | Reducing Agent | 48. (D12) | 0.001 M (1 mM) |

Additive Screen contains ninety-six unique reagents beginning at position A1.
To determine the formulation of each reagent, simply read across the page.

Lot 2428** where ** ≥ 18

| Tube # | Salt | Tube # | Classification | Tube # | Suggested Drop Concentration |
|-----------|---|-----------|-----------------------|-----------|------------------------------|
| 49. (E1) | 0.1 M Ethylenediaminetetraacetic acid disodium salt dihydrate | 49. (E1) | Chelating Agent | 49. (E1) | 0.01 M (10 mM) |
| 50. (E2) | 5% w/v Polyvinylpyrrolidone K15 | 50. (E2) | Polymer | 50. (E2) | 0.5% |
| 51. (E3) | 30% w/v Dextran sulfate sodium salt (M _r 5,000) | 51. (E3) | Polymer | 51. (E3) | 3.0% |
| 52. (E4) | 40% v/v Pentaerythritol ethoxylate (3/4 EO/OH) | 52. (E4) | Polymer | 52. (E4) | 4.0% |
| 53. (E5) | 10% w/v Polyethylene glycol 3,350 | 53. (E5) | Polymer | 53. (E5) | 1.0% |
| 54. (E6) | 30% w/v D-(+)-Glucose monohydrate | 54. (E6) | Carbohydrate | 54. (E6) | 3.0% |
| 55. (E7) | 30% w/v Sucrose | 55. (E7) | Carbohydrate | 55. (E7) | 3.0% |
| 56. (E8) | 30% w/v Xylitol | 56. (E8) | Carbohydrate | 56. (E8) | 3.0% |
| 57. (E9) | 30% w/v D-Sorbitol | 57. (E9) | Carbohydrate | 57. (E9) | 3.0% |
| 58. (E10) | 12% w/v myo-Inositol | 58. (E10) | Carbohydrate | 58. (E10) | 1.2% |
| 59. (E11) | 30% w/v D-(+)-Trehalose dihydrate | 59. (E11) | Carbohydrate | 59. (E11) | 3.0% |
| 60. (E12) | 30% w/v D-(+)-Galactose | 60. (E12) | Carbohydrate | 60. (E12) | 3.0% |
| 61. (F1) | 30% v/v Ethylene glycol | 61. (F1) | Polyol | 61. (F1) | 3.0% |
| 62. (F2) | 30% v/v Glycerol | 62. (F2) | Polyol | 62. (F2) | 3.0% |
| 63. (F3) | 3.0 M NDSB-195 | 63. (F3) | Non-detergent | 63. (F3) | 0.3 M (300 mM) |
| 64. (F4) | 2.0 M NDSB-201 | 64. (F4) | Non-detergent | 64. (F4) | 0.2 M (200 mM) |
| 65. (F5) | 2.0 M NDSB-211 | 65. (F5) | Non-detergent | 65. (F5) | 0.2 M (200 mM) |
| 66. (F6) | 2.0 M NDSB-221 | 66. (F6) | Non-detergent | 66. (F6) | 0.2 M (200 mM) |
| 67. (F7) | 1.0 M NDSB-256 | 67. (F7) | Non-detergent | 67. (F7) | 0.1 M (200 mM) |
| 68. (F8) | 0.5% w/v 1,2,3-Heptanetriol | 68. (F8) | Amphiphile | 68. (F8) | 0.05% |
| 69. (F9) | 20% w/v Benzamidine hydrochloride | 69. (F9) | Amphiphile | 69. (F9) | 2.0% |
| 70. (F10) | 5% w/v n-dodecyl-N,N-dimethylamine-N-oxide, (LDAO, DDAO) | 70. (F10) | Detergent | 70. (F10) | 0.5% |
| 71. (F11) | 5% w/v n-Octyl-β-D-glucoside | 71. (F11) | Detergent | 71. (F11) | 0.5% |
| 72. (F12) | 5% w/v n-Dodecyl-β-D-maltoside | 72. (F12) | Detergent | 72. (F12) | 0.5% |
| 73. (G1) | 30% w/v Trimethylamine N-oxide dihydrate | 73. (G1) | Osmolyte | 73. (G1) | 3.0% |
| 74. (G2) | 30% w/v 1,6-Hexanediol | 74. (G2) | Organic, Non-volatile | 74. (G2) | 3.0% |
| 75. (G3) | 30% v/v (+/-)-2-Methyl-2,4-pentanediol | 75. (G3) | Organic, Non-volatile | 75. (G3) | 3.0% |
| 76. (G4) | 50% v/v Polyethylene glycol 400 | 76. (G4) | Organic, Non-volatile | 76. (G4) | 5.0% |
| 77. (G5) | 50% v/v Jeffamine® M-600® pH 7.0 | 77. (G5) | Organic, Non-volatile | 77. (G5) | 5.0% |
| 78. (G6) | 40% v/v 2,5-Hexanediol (mixture of isomers) | 78. (G6) | Organic, Non-volatile | 78. (G6) | 4.0% |
| 79. (G7) | 40% v/v (±)-1,3-Butanediol | 79. (G7) | Organic, Non-volatile | 79. (G7) | 4.0% |
| 80. (G8) | 40% v/v Polypropylene glycol P 400 | 80. (G8) | Organic, Non-volatile | 80. (G8) | 4.0% |
| 81. (G9) | 30% v/v 1,4-Dioxane | 81. (G9) | Organic, Volatile | 81. (G9) | 3.0% |
| 82. (G10) | 30% v/v Ethanol | 82. (G10) | Organic, Volatile | 82. (G10) | 3.0% |
| 83. (G11) | 30% v/v 2-Propanol | 83. (G11) | Organic, Volatile | 83. (G11) | 3.0% |
| 84. (G12) | 30% v/v Methanol | 84. (G12) | Organic, Volatile | 84. (G12) | 3.0% |
| 85. (H1) | 10% v/v 1,2-Butanediol | 85. (H1) | Organic, Volatile | 85. (H1) | 1.0% |
| 86. (H2) | 40% v/v tert-Butanol | 86. (H2) | Organic, Volatile | 86. (H2) | 4.0% |
| 87. (H3) | 40% v/v 1,3-Propanediol | 87. (H3) | Organic, Volatile | 87. (H3) | 4.0% |
| 88. (H4) | 40% v/v Acetonitrile | 88. (H4) | Organic, Volatile | 88. (H4) | 4.0% |
| 89. (H5) | 40% v/v Formamide | 89. (H5) | Organic, Volatile | 89. (H5) | 4.0% |
| 90. (H6) | 40% v/v 1-Propanol | 90. (H6) | Organic, Volatile | 90. (H6) | 4.0% |
| 91. (H7) | 5% v/v Ethyl acetate | 91. (H7) | Organic, Volatile | 91. (H7) | 0.5% |
| 92. (H8) | 40% v/v Acetone | 92. (H8) | Organic, Volatile | 92. (H8) | 4.0% |
| 93. (H9) | 0.25% v/v Dichloromethane | 93. (H9) | Organic, Volatile | 93. (H9) | 0.025% |
| 94. (H10) | 7% v/v 1-Butanol | 94. (H10) | Organic, Volatile | 94. (H10) | 0.7% |
| 95. (H11) | 40% v/v 2,2,2-Trifluoroethanol | 95. (H11) | Organic, Volatile | 95. (H11) | 4.0% |
| 96. (H12) | 40% v/v 1,1,1,3,3,3-Hexafluoro-2-propanol | 96. (H12) | Organic, Volatile | 96. (H12) | 4.0% |

Additive Screen contains ninety-six unique reagents beginning at position A1.
To determine the formulation of each reagent, simply read across the page.

Lot 2428** where ** ≥18

Sample: _____ Sample Concentration: _____
 Sample Buffer: _____ Date: _____
 Reservoir Volume: _____ Temperature: _____
 Drop Volume: Total _____ μ l Sample _____ μ l Reservoir _____ μ l Additive _____ μ l

1 Clear Drop
 2 Phase Separation
 3 Regular Granular Precipitate
 4 Birefringent Precipitate or Microcrystals

5 Posettes or Spherulites
 6 Needles (1D Growth)
 7 Plates (2D Growth)
 8 Single Crystals (3D Growth < 0.2 mm)
 9 Single Crystals (3D Growth > 0.2 mm)

Additive Screen™ - HR2-428 Scoring Sheet

| | | Date: | Date: | Date: | Date: |
|-----------|---|--------------------|-------|-------|-------|
| 1. (A1) | 0.1 M Barium chloride dihydrate | Multivalent | | | |
| 2. (A2) | 0.1 M Cadmium chloride hydrate | Multivalent | | | |
| 3. (A3) | 0.1 M Calcium chloride dihydrate | Multivalent | | | |
| 4. (A4) | 0.1 M Cobalt(II) chloride hexahydrate | Multivalent | | | |
| 5. (A5) | 0.1 M Copper(II) chloride dihydrate | Multivalent | | | |
| 6. (A6) | 0.1 M Magnesium chloride hexahydrate | Multivalent | | | |
| 7. (A7) | 0.1 M Manganese(II) chloride tetrahydrate | Multivalent | | | |
| 8. (A8) | 0.1 M Strontium chloride hexahydrate | Multivalent | | | |
| 9. (A9) | 0.1 M Yttrium(III) chloride hexahydrate | Multivalent | | | |
| 10. (A10) | 0.1 M Zinc chloride | Multivalent | | | |
| 11. (A11) | 0.1 M Iron(III) chloride hexahydrate | Multivalent | | | |
| 12. (A12) | 0.1 M Nickel(II) chloride hexahydrate | Multivalent | | | |
| 13. (B1) | 0.1 M Chromium(III) chloride hexahydrate | Multivalent | | | |
| 14. (B2) | 0.1 M Praseodymium(III) acetate hydrate | Multivalent | | | |
| 15. (B3) | 1.0 M Ammonium sulfate | Salt | | | |
| 16. (B4) | 1.0 M Potassium chloride | Salt | | | |
| 17. (B5) | 1.0 M Lithium chloride | Salt | | | |
| 18. (B6) | 2.0 M Sodium chloride | Salt | | | |
| 19. (B7) | 0.5 M Sodium fluoride | Salt | | | |
| 20. (B8) | 1.0 M Sodium iodide | Salt | | | |
| 21. (B9) | 2.0 M Sodium thiocyanate | Salt | | | |
| 22. (B10) | 1.0 M Potassium sodium tartrate tetrahydrate | Salt | | | |
| 23. (B11) | 1.0 M Sodium citrate tribasic dihydrate | Salt | | | |
| 24. (B12) | 1.0 M Cesium chloride | Salt | | | |
| 25. (C1) | 1.0 M Sodium malonate pH 7.0 | Salt | | | |
| 26. (C2) | 0.1 M L-Proline | Amino Acid | | | |
| 27. (C3) | 0.1 M Phenol | Dissociating Agent | | | |
| 28. (C4) | 30% v/v Dimethyl sulfoxide | Dissociating Agent | | | |
| 29. (C5) | 0.1 M Sodium bromide | Dissociating Agent | | | |
| 30. (C6) | 30% w/v 6-Aminohexanoic acid | Linker | | | |
| 31. (C7) | 30% w/v 1,5-Diaminopentane dihydrochloride | Linker | | | |
| 32. (C8) | 30% w/v 1,6-Diaminohexane | Linker | | | |
| 33. (C9) | 30% w/v 1,8-Diaminooctane | Linker | | | |
| 34. (C10) | 1.0 M Glycine | Linker | | | |
| 35. (C11) | 0.3 M Glycyl-glycyl-glycine | Linker | | | |
| 36. (C12) | 0.1 M Taurine | Linker | | | |
| 37. (D1) | 0.1 M Betaine hydrochloride | Linker | | | |
| 38. (D2) | 0.1 M Spermidine | Polyamine | | | |
| 39. (D3) | 0.1 M Spermine tetrahydrochloride | Polyamine | | | |
| 40. (D4) | 0.1 M Hexamine cobalt(III) chloride | Polyamine | | | |
| 41. (D5) | 0.1 M Sarcosine | Polyamine/Osmolyte | | | |
| 42. (D6) | 0.1 M Trimethylamine hydrochloride | Chaotrope | | | |
| 43. (D7) | 1.0 M Guanidine hydrochloride | Chaotrope | | | |
| 44. (D8) | 0.1 M Urea | Chaotrope | | | |
| 45. (D9) | 0.1 M β -Nicotinamide adenine dinucleotide hydrate | Co-factor | | | |
| 46. (D10) | 0.1 M Adenosine-5'-triphosphate disodium salt hydrate | Co-factor | | | |
| 47. (D11) | 0.1 M TCEP hydrochloride | Reducing agent | | | |
| 48. (D12) | 0.01 M GSH (L-Glutathione reduced), 0.01 M GSSG (L-Glutathione oxidized) | Reducing agent | | | |



Solutions for Crystal Growth

34 Journey
 Aliso Viejo, CA 92656-3317 U.S.A.
 Tel: (949) 425-1321 • Fax: (949) 425-1611
 e-mail: tech@hrcmail.com
 Website: www.hamptonresearch.com

Sample: _____ Sample Concentration: _____
 Sample Buffer: _____ Date: _____
 Reservoir Volume: _____ Temperature: _____
 Drop Volume: Total _____ μ l Sample _____ μ l Reservoir _____ μ l Additive _____ μ l

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals
- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2 mm)
- 9 Single Crystals (3D Growth > 0.2 mm)

Additive Screen™ - HR2-428 Scoring Sheet

| | | Date: | Date: | Date: | Date: |
|-----------|---|-----------------------|-------|-------|-------|
| 49. (E1) | 0.1 M Ethylenediaminetetraacetic acid disodium salt dihydrate | Chelating Agent | | | |
| 50. (E2) | 5% w/v Polyvinylpyrrolidone K15 | Polymer | | | |
| 51. (E3) | 30% w/v Dextran sulfate sodium salt (Mr 5,000) | Polymer | | | |
| 52. (E4) | 40% v/v Pentaerythritol ethoxylate (3/4 EO/OH) | Polymer | | | |
| 53. (E5) | 10% w/v Polyethylene glycol 3,350 | Polymer | | | |
| 54. (E6) | 30% w/v D-(+)-Glucose monohydrate | Carbohydrate | | | |
| 55. (E7) | 30% w/v Sucrose | Carbohydrate | | | |
| 56. (E8) | 30% w/v Xylitol | Carbohydrate | | | |
| 57. (E9) | 30% w/v D-Sorbitol | Carbohydrate | | | |
| 58. (E10) | 12% w/v myo-Inositol | Carbohydrate | | | |
| 59. (E11) | 30% w/v D-(+)-Trehalose dihydrate | Carbohydrate | | | |
| 60. (E12) | 30% w/v D-(+)-Galactose | Carbohydrate | | | |
| 61. (F1) | 30% v/v Ethylene glycol | Polyol | | | |
| 62. (F2) | 30% v/v Glycerol | Polyol | | | |
| 63. (F3) | 3.0 M NDSB-195 | Non-detergent | | | |
| 64. (F4) | 2.0 M NDSB-201 | Non-detergent | | | |
| 65. (F5) | 2.0 M NDSB-211 | Non-detergent | | | |
| 66. (F6) | 2.0 M NDSB-221 | Non-detergent | | | |
| 67. (F7) | 1.0 M NDSB-256 | Non-detergent | | | |
| 68. (F8) | 0.5% w/v 1,2,3-Heptanetriol | Amphiphile | | | |
| 69. (F9) | 20% w/v Benzamidine hydrochloride | Amphiphile | | | |
| 70. (F10) | 5% w/v n-dodecyl-N,N-dimethylamine-N-oxide, (LDAO, DDAO) | Detergent | | | |
| 71. (F11) | 5% w/v n-Octyl- β -D-glucoside | Detergent | | | |
| 72. (F12) | 5% w/v n-Dodecyl- β -D-maltoside | Detergent | | | |
| 73. (G1) | 30% w/v Trimethylamine N-oxide dihydrate | Osmolyte | | | |
| 74. (G2) | 30% w/v 1,6-Hexanediol | Organic, Non-volatile | | | |
| 75. (G3) | 30% v/v (+/-)-2-Methyl-2,4-pentanediol | Organic, Non-volatile | | | |
| 76. (G4) | 50% v/v Polyethylene glycol 400 | Organic, Non-volatile | | | |
| 77. (G5) | 50% v/v Jeffamine® M-600® pH 7.0 | Organic, Non-volatile | | | |
| 78. (G6) | 40% v/v 2,5-Hexanediol | Organic, Non-volatile | | | |
| 79. (G7) | 40% v/v (\pm)-1,3-Butanediol | Organic, Non-volatile | | | |
| 80. (G8) | 40% v/v Polypropylene glycol P 400 | Organic, Non-volatile | | | |
| 81. (G9) | 30% v/v 1,4-Dioxane | Organic, Volatile | | | |
| 82. (G10) | 30% v/v Ethanol | Organic, Volatile | | | |
| 83. (G11) | 30% v/v 2-Propanol | Organic, Volatile | | | |
| 84. (G12) | 30% v/v Methanol | Organic, Volatile | | | |
| 85. (H1) | 10% v/v 1,2-Butanediol | Organic, Volatile | | | |
| 86. (H2) | 40% v/v tert-Butanol | Organic, Volatile | | | |
| 87. (H3) | 40% v/v 1,3-Propanediol | Organic, Volatile | | | |
| 88. (H4) | 40% v/v Acetonitrile | Organic, Volatile | | | |
| 89. (H5) | 40% v/v Formamide | Organic, Volatile | | | |
| 90. (H6) | 40% v/v 1-Propanol | Organic, Volatile | | | |
| 91. (H7) | 5% v/v Ethyl acetate | Organic, Volatile | | | |
| 92. (H8) | 40% v/v Acetone | Organic, Volatile | | | |
| 93. (H9) | 0.25% v/v Dichloromethane | Organic, Volatile | | | |
| 94. (H10) | 7% v/v 1-Butanol | Organic, Volatile | | | |
| 95. (H11) | 40% v/v 2,2,2-Trifluoroethanol | Organic, Volatile | | | |
| 96. (H12) | 40% v/v 1,1,1,3,3,3-Hexafluoro-2-propanol | Organic, Volatile | | | |



Solutions for Crystal Growth

34 Journey
 Also Viejo, CA 92586-3317 U.S.A.
 Tel: (949) 425-1321 • Fax: (949) 425-1611
 e-mail: tech@hrcr.com
 Website: www.hamptonresearch.com

© 1991-2019 Hampton Research Corp. all rights reserved
 Printed in the United States of America. This guide or parts thereof may not
 be reproduced in any form without the written permission of the publisher.