

# Frequently Asked Questions

HAMPTON  
RESEARCH

*Solutions for Crystal Growth*

## Additive Screen and Additive Screen HT

## HR2-138 & HR2-428

### **When and why did reagent 89(H5) change from 40% v/v $\gamma$ -Butyrolactone to 40% v/v Formamide?**

$\gamma$ -Butyrolactone became a United States of America Drug Enforcement Agency (DEA) controlled substance, making it difficult or impossible to obtain for many labs, since a DEA permit is required to order and use this chemical. Hampton Research substituted Formamide for reagent 89(H5) during June 2008.

Additive Screen kits (HR2-428) with the lot number 242801, 242802, 242803, 242804, 242805, 242806, 242807, and 242808 contain 40% v/v  $\gamma$ -Butyrolactone as reagent 89(H5). Additive Screen kits (HR2-428) with lot number 242809 to 2428\*\* where \*\*  $\geq 09$  contain 40% v/v Formamide as reagent 89(H5).

Additive Screen HT kits (HR2-138) with the lot number 213801, 213802, and 213803 contain 40% v/v  $\gamma$ -Butyrolactone as reagent 89(H5). Additive Screen HT kits (HR2-138) with the lot number 213804 to 2138\*\* were \*\*  $\geq 4$  contain 40% v/v Formamide as reagent 89(H5).

### **When and why did reagent 85(H1) change from 40% v/v 1,4-Butanediol replaced with 10% v/v 1,2-Butanediol.**

1,4-Butanediol became a United States of America Drug Enforcement Agency (DEA) controlled substance, making it difficult or impossible to obtain for many labs, since a DEA permit is required to order and use this chemical. Hampton Research substituted 1,2-Butanediol for reagent 85(H1) during April 2010.

Additive Screen kits (HR2-428) with the lot number 242801, 242802, 242803, 242804, 242805, 242806, 242807, 242808, 242809 and 242810 contain 40% v/v 1,4-Butanediol as reagent 85(H1). Additive Screen kits (HR2-428) with lot number 242811 to 2428\*\* where \*\*  $\geq 11$  contain 10% v/v 1,2-Butanediol as reagent 85(H1).

Additive Screen HT kits (HR2-138) with the lot number 213801, 213802, 213803, 218304, and 218305 contain 40% v/v 1,4-Butanediol as reagent 85(H1). Additive Screen HT kits (HR2-138) with the lot number 213806 to 2138\*\* were \*\*  $\geq 6$  contain 10% v/v 1,2-Butanediol for reagent 85(H1).

### **Is the 2,5-hexanediol a racemic mixture (mixture of L and R isomers), or the pure enantiomer (L or R form)?**

The 2,5-hexanediol used in the Additive Screen and Additive Screen HT is a racemic mixture (mixture of L and R isomers).

Hampton Research  
34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

Technical Support e-mail: tech@hrmail.com

Website: www.hamptonresearch.com