

Well #	Salt	Well #	Buffer ◇	Well #	Precipitant
1. (A1)	0.01 M Magnesium chloride hexahydrate	1. (A1)	0.05 M MES monohydrate pH 5.6	1. (A1)	1.8 M Lithium sulfate monohydrate
2. (A2)	0.01 M Magnesium acetate tetrahydrate	2. (A2)	0.05 M MES monohydrate pH 5.6	2. (A2)	2.5 M Ammonium sulfate
3. (A3)	0.1 M Magnesium acetate tetrahydrate	3. (A3)	0.05 M MES monohydrate pH 5.6	3. (A3)	20% v/v (+/-)-2-Methyl-2,4-pentanediol
4. (A4)	0.2 M Potassium chloride, 0.01 M Magnesium sulfate heptahydrate	4. (A4)	0.05 M MES monohydrate pH 5.6	4. (A4)	10% v/v Polyethylene glycol 400
5. (A5)	0.2 M Potassium chloride, 0.01 M Magnesium chloride hexahydrate	5. (A5)	0.05 M MES monohydrate pH 5.6	5. (A5)	5% w/v Polyethylene glycol 8,000
6. (A6)	0.1 M Ammonium sulfate, 0.01 M Magnesium chloride hexahydrate	6. (A6)	0.05 M MES monohydrate pH 5.6	6. (A6)	20% w/v Polyethylene glycol 8,000
7. (A7)	0.02 M Magnesium chloride hexahydrate	7. (A7)	0.05 M MES monohydrate pH 6.0	7. (A7)	15% v/v 2-Propanol
8. (A8)	0.1 M Ammonium acetate, 0.005 M Magnesium sulfate heptahydrate	8. (A8)	0.05 M MES monohydrate pH 6.0	8. (A8)	0.6 M Sodium chloride
9. (A9)	0.1 M Potassium chloride, 0.01 M Magnesium chloride hexahydrate	9. (A9)	0.05 M MES monohydrate pH 6.0	9. (A9)	10% v/v Polyethylene glycol 400
10. (A10)	0.005 M Magnesium sulfate heptahydrate	10. (A10)	0.05 M MES monohydrate pH 6.0	10. (A10)	5% w/v Polyethylene glycol 4,000
11. (A11)	0.01 M Magnesium chloride hexahydrate	11. (A11)	0.05 M Sodium cacodylate trihydrate pH 6.0	11. (A11)	1.0 M Lithium sulfate monohydrate
12. (A12)	0.01 M Magnesium sulfate heptahydrate	12. (A12)	0.05 M Sodium cacodylate trihydrate pH 6.0	12. (A12)	1.8 M Lithium sulfate monohydrate
13. (B1)	0.015 M Magnesium acetate tetrahydrate	13. (B1)	0.05 M Sodium cacodylate trihydrate pH 6.0	13. (B1)	1.7 M Ammonium sulfate
14. (B2)	0.1 M Potassium chloride, 0.025 M Magnesium chloride hexahydrate	14. (B2)	0.05 M Sodium cacodylate trihydrate pH 6.0	14. (B2)	15% v/v 2-Propanol
15. (B3)	0.04 M Magnesium chloride hexahydrate	15. (B3)	0.05 M Sodium cacodylate trihydrate pH 6.0	15. (B3)	5% v/v (+/-)-2-Methyl-2,4-pentanediol
16. (B4)	0.04 M Magnesium sulfate tetrahydrate	16. (B4)	0.05 M Sodium cacodylate trihydrate pH 6.0	16. (B4)	30% v/v (+/-)-2-Methyl-2,4-pentanediol
17. (B5)	0.2 M Potassium chloride, 0.01 M Calcium chloride dihydrate	17. (B5)	0.05 M Sodium cacodylate trihydrate pH 6.0	17. (B5)	10% w/v Polyethylene glycol 4,000
18. (B6)	0.01 M Magnesium acetate tetrahydrate	18. (B6)	0.05 M Sodium cacodylate trihydrate pH 6.5	18. (B6)	1.3 M Lithium sulfate monohydrate
19. (B7)	0.01 M Magnesium sulfate heptahydrate	19. (B7)	0.05 M Sodium cacodylate trihydrate pH 6.5	19. (B7)	2.0 M Ammonium sulfate
20. (B8)	0.1 M Ammonium acetate, 0.015 M Magnesium acetate tetrahydrate	20. (B8)	0.05 M Sodium cacodylate trihydrate pH 6.5	20. (B8)	10% v/v 2-Propanol
21. (B9)	0.2 M Potassium chloride, 0.005 M Magnesium chloride hexahydrate	21. (B9)	0.05 M Sodium cacodylate trihydrate pH 6.5	21. (B9)	0.9 M 1,6-Hexanediol
22. (B10)	0.08 M Magnesium acetate tetrahydrate	22. (B10)	0.05 M Sodium cacodylate trihydrate pH 6.5	22. (B10)	15% v/v Polyethylene glycol 400
23. (B11)	0.2 M Potassium chloride, 0.01 Magnesium chloride hexahydrate	23. (B11)	0.05 M Sodium cacodylate trihydrate pH 6.5	23. (B11)	10% w/v Polyethylene glycol 4,000
24. (B12)	0.2 M Ammonium acetate, 0.01 M Calcium chloride dihydrate	24. (B12)	0.05 M Sodium cacodylate trihydrate pH 6.5	24. (B12)	10% w/v Polyethylene glycol 4,000
25. (C1)	0.08 M Magnesium acetate tetrahydrate	25. (C1)	0.05 M Sodium cacodylate trihydrate pH 6.5	25. (C1)	30% w/v Polyethylene glycol 4,000
26. (C2)	0.2 M Potassium chloride, 0.1 M Magnesium acetate tetrahydrate	26. (C2)	0.05 M Sodium cacodylate trihydrate pH 6.5	26. (C2)	10% w/v Polyethylene glycol 8,000
27. (C3)	0.2 M Ammonium acetate, 0.01 M Magnesium acetate tetrahydrate	27. (C3)	0.05 M Sodium cacodylate trihydrate pH 6.5	27. (C3)	30% v/v Polyethylene glycol 8,000
28. (C4)	0.05 M Magnesium sulfate hydrate	28. (C4)	0.05 M HEPES sodium pH 7.0	28. (C4)	1.6 M Lithium sulfate monohydrate
29. (C5)	0.01 M Magnesium chloride hexahydrate	29. (C5)	0.05 M HEPES sodium pH 7.0	29. (C5)	4.0 M Lithium chloride
30. (C6)	0.01 M Magnesium chloride hexahydrate	30. (C6)	0.05 M HEPES sodium pH 7.0	30. (C6)	1.6 M Ammonium sulfate
31. (C7)	0.005 M Magnesium chloride hexahydrate	31. (C7)	0.05 M HEPES sodium pH 7.0	31. (C7)	25% v/v Polyethylene glycol monomethyl ether 550
32. (C8)	0.2 M Potassium chloride, 0.01 M Magnesium chloride hexahydrate	32. (C8)	0.05 M HEPES sodium pH 7.0	32. (C8)	1.7 M 1,6-Hexanediol
33. (C9)	0.2 M Ammonium chloride, 0.01 M Magnesium chloride hexahydrate	33. (C9)	0.05 M HEPES sodium pH 7.0	33. (C9)	2.5 M 1,6-Hexanediol
34. (C10)	0.1 M Potassium chloride, 0.005 M Magnesium sulfate hydrate	34. (C10)	0.05 M HEPES sodium pH 7.0	34. (C10)	15% v/v (+/-)-2-Methyl-2,4-pentanediol
35. (C11)	0.1 M Potassium chloride, 0.01 M Magnesium chloride hexahydrate	35. (C11)	0.05 M HEPES sodium pH 7.0	35. (C11)	5% v/v Polyethylene glycol 400
36. (C12)	0.1 M Potassium chloride, 0.01 M Calcium chloride dihydrate	36. (C12)	0.05 M HEPES sodium pH 7.0	36. (C12)	10% v/v Polyethylene glycol 400
37. (D1)	0.2 M Potassium chloride, 0.025 M Magnesium sulfate hydrate	37. (D1)	0.05 M HEPES sodium pH 7.0	37. (D1)	20% v/v Polyethylene glycol 200
38. (D2)	0.2 M Ammonium acetate, 0.15 M Magnesium acetate tetrahydrate	38. (D2)	0.05 M HEPES sodium pH 7.0	38. (D2)	5% w/v Polyethylene glycol 4,000
39. (D3)	0.1 M Ammonium acetate, 0.02 M Magnesium chloride hexahydrate	39. (D3)	0.05 M HEPES sodium pH 7.0	39. (D3)	5% w/v Polyethylene glycol 8,000
40. (D4)	0.01 M Magnesium chloride hexahydrate	40. (D4)	0.05 M TRIS hydrochloride pH 7.5	40. (D4)	1.6 M Ammonium sulfate
41. (D5)	0.1 M Potassium chloride, 0.015 M Magnesium chloride hexahydrate	41. (D5)	0.05 M TRIS hydrochloride pH 7.5	41. (D5)	10% v/v Polyethylene glycol monomethyl ether 550
42. (D6)	0.01 M Magnesium chloride hexahydrate	42. (D6)	0.05 M TRIS hydrochloride pH 7.5	42. (D6)	5% v/v 2-Propanol
43. (D7)	0.05 M Ammonium acetate, 0.01 M Magnesium chloride hexahydrate	43. (D7)	0.05 M TRIS hydrochloride pH 7.5	43. (D7)	10% v/v (+/-)-2-Methyl-2,4-pentanediol
44. (D8)	0.2 M Potassium chloride, 0.05 M Magnesium chloride hexahydrate	44. (D8)	0.05 M TRIS hydrochloride pH 7.5	44. (D8)	10% w/v Polyethylene glycol 4,000
45. (D9)	0.025 M Magnesium sulfate hydrate	45. (D9)	0.05 M TRIS hydrochloride pH 8.5	45. (D9)	1.8 M Ammonium sulfate
46. (D10)	0.005 M Magnesium sulfate hydrate	46. (D10)	0.05 M TRIS hydrochloride pH 8.5	46. (D10)	2.9 M 1,6-Hexanediol
47. (D11)	0.1 M Potassium chloride, 0.01 M Magnesium chloride hexahydrate	47. (D11)	0.05 M TRIS hydrochloride pH 8.5	47. (D11)	30% v/v Polyethylene glycol 400
48. (D12)	0.2 M Ammonium chloride, 0.01 M Calcium chloride dihydrate	48. (D12)	0.05 M TRIS hydrochloride pH 8.5	48. (D12)	30% w/v Polyethylene glycol 4,000

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

Natrix HT contains ninety-six unique reagents. To determine the formulation of each reagent, simply read across the page.



Solutions for Crystal Growth

Well #	Salt	Well #	Buffer ◊	Well #	Precipitant	Well #	Additive
49.(E1)	0.04 M Lithium chloride, 0.02 M Magnesium chloride hexahydrate	49.(E1)	0.04 M Sodium cacodylate trihydrate pH 5.5	49.(E1)	30% v/v (+/-)-2-Methyl-2,4-pentanediol	49.(E1)	0.002 M Hexamine cobalt(III) chloride
50.(E2)	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate	50.(E2)	0.04 M Sodium cacodylate trihydrate pH 5.5	50.(E2)	35% v/v (+/-)-2-Methyl-2,4-pentanediol	50.(E2)	0.002 M Hexamine cobalt(III) chloride
51.(E3)	0.012 M Sodium chloride, 0.08 M Potassium chloride	51.(E3)	0.04 M Sodium cacodylate trihydrate pH 5.5	51.(E3)	45% v/v (+/-)-2-Methyl-2,4-pentanediol	51.(E3)	0.002 M Hexamine cobalt(III) chloride
52.(E4)	0.02 M Magnesium chloride hexahydrate	52.(E4)	0.04 M Sodium cacodylate trihydrate pH 5.5	52.(E4)	40% v/v (+/-)-2-Methyl-2,4-pentanediol	52.(E4)	0.002 M Hexamine cobalt(III) chloride
53.(E5)	0.002 M Calcium chloride dihydrate	53.(E5)	0.05 M Sodium cacodylate trihydrate pH 6.0	53.(E5)	1.8 M Ammonium sulfate	53.(E5)	0.0005 M Spermine
54.(E6)		54.(E6)	0.05 M Sodium cacodylate trihydrate pH 6.0	54.(E6)	35% v/v Tacsimate™ pH 6.0	54.(E6)	0.001 M Spermine
55.(E7)	0.1 M Sodium chloride	55.(E7)	0.05 M Sodium cacodylate trihydrate pH 6.0	55.(E7)	10% w/v Polyethylene glycol 4,000	55.(E7)	0.0005 M Spermine
56.(E8)	0.05 M Potassium chloride	56.(E8)	0.05 M Sodium cacodylate trihydrate pH 6.0	56.(E8)	10% w/v Polyethylene glycol 8,000	56.(E8)	0.0005 M Spermine, 0.0005 M L-Argininamide dihydrochloride
57.(E9)	0.1 M Potassium chloride	57.(E9)	0.05 M Sodium cacodylate trihydrate pH 6.0	57.(E9)	16% w/v Polyethylene glycol 1,000	57.(E9)	0.0005 M Spermine
58.(E10)	0.005 M Magnesium chloride hexahydrate, 0.002 M Calcium chloride dihydrate	58.(E10)	0.05 M Sodium cacodylate trihydrate pH 6.0	58.(E10)	15% v/v 2-Propanol	58.(E10)	0.001 M Spermine
59.(E11)	0.075 M Sodium chloride, 0.002 M Calcium chloride dihydrate	59.(E11)	0.05 M Sodium cacodylate trihydrate pH 6.0	59.(E11)	30% w/v 1,6-Hexanediol	59.(E11)	0.0005 M Spermine
60.(E12)	0.02 M Magnesium sulfate hydrate, 0.002 M Cobalt(II) chloride hexahydrate	60.(E12)	0.05 M Sodium cacodylate trihydrate pH 6.0	60.(E12)	25% v/v (+/-)-2-Methyl-2,4-pentanediol	60.(E12)	0.0005 M Spermine
61.(F1)		61.(F1)	0.05 M Sodium cacodylate trihydrate pH 6.0	61.(F1)	30% v/v (+/-)-2-Methyl-2,4-pentanediol	61.(F1)	
62.(F2)	0.08 M Sodium chloride, 0.012 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	62.(F2)	0.04 M Sodium cacodylate trihydrate pH 6.0	62.(F2)	30% v/v (+/-)-2-Methyl-2,4-pentanediol	62.(F2)	0.012 M Spermine tetrahydrochloride
63.(F3)	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate	63.(F3)	0.04 M Sodium cacodylate trihydrate pH 6.0	63.(F3)	35% v/v (+/-)-2-Methyl-2,4-pentanediol	63.(F3)	0.012 M Spermine tetrahydrochloride
64.(F4)	0.08 M Strontium chloride hexahydrate	64.(F4)	0.04 M Sodium cacodylate trihydrate pH 6.0	64.(F4)	35% v/v (+/-)-2-Methyl-2,4-pentanediol	64.(F4)	0.012 M Spermine tetrahydrochloride
65.(F5)	0.08 M Potassium chloride, 0.02 M Barium chloride dihydrate	65.(F5)	0.04 M Sodium cacodylate trihydrate pH 6.0	65.(F5)	40% v/v (+/-)-2-Methyl-2,4-pentanediol	65.(F5)	0.012 M Spermine tetrahydrochloride
66.(F6)	0.08 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	66.(F6)	0.04 M Sodium cacodylate trihydrate pH 6.0	66.(F6)	45% v/v (+/-)-2-Methyl-2,4-pentanediol	66.(F6)	0.012 M Spermine tetrahydrochloride
67.(F7)	0.08 M Sodium chloride	67.(F7)	0.04 M Sodium cacodylate trihydrate pH 6.0	67.(F7)	45% v/v (+/-)-2-Methyl-2,4-pentanediol	67.(F7)	0.012 M Spermine tetrahydrochloride
68.(F8)	0.08 M Sodium chloride, 0.02 M Barium chloride dihydrate	68.(F8)	0.04 M Sodium cacodylate trihydrate pH 6.0	68.(F8)	45% v/v (+/-)-2-Methyl-2,4-pentanediol	68.(F8)	0.012 M Spermine tetrahydrochloride
69.(F9)	0.012 M Sodium chloride, 0.08 M Potassium chloride	69.(F9)	0.04 M Sodium cacodylate trihydrate pH 6.0	69.(F9)	50% v/v (+/-)-2-Methyl-2,4-pentanediol	69.(F9)	0.012 M Spermine tetrahydrochloride
70.(F10)	0.08 M Potassium chloride	70.(F10)	0.04 M Sodium cacodylate trihydrate pH 6.0	70.(F10)	55% v/v (+/-)-2-Methyl-2,4-pentanediol	70.(F10)	0.012 M Spermine tetrahydrochloride
71.(F11)	0.018 M Magnesium chloride hexahydrate	71.(F11)	0.05 M Sodium cacodylate trihydrate pH 6.5	71.(F11)	10% v/v 2-Propanol	71.(F11)	0.003 M Spermine
72.(F12)	0.02 M Magnesium chloride hexahydrate	72.(F12)	0.05 M MOPS pH 7.0	72.(F12)	2.0 M Ammonium sulfate	72.(F12)	0.0005 M Spermine
73.(G1)		73.(G1)	0.05 M HEPES sodium pH 7.0	73.(G1)	40% v/v Tacsimate™ pH 7.0	73.(G1)	0.002 M Spermine, 0.002 M Hexamine cobalt(III) chloride
74.(G2)	0.02 M Magnesium chloride hexahydrate	74.(G2)	0.05 M MOPS pH 7.0	74.(G2)	55% v/v Tacsimate™ pH 7.0	74.(G2)	0.002 M Hexamine cobalt(III) chloride
75.(G3)	0.02 M Magnesium chloride hexahydrate	75.(G3)	0.05 M Sodium cacodylate trihydrate pH 7.0	75.(G3)	15% v/v 2-Propanol	75.(G3)	0.001 M Hexamine cobalt(III) chloride, 0.001 M Spermine
76.(G4)	0.005 M Magnesium chloride hexahydrate	76.(G4)	0.05 M MOPS pH 7.0	76.(G4)	25% v/v 1,4-Dioxane	76.(G4)	0.001 M Spermine
77.(G5)	0.01 M Magnesium chloride hexahydrate, 0.002 M Barium chloride dihydrate	77.(G5)	0.05 M MOPS pH 7.0	77.(G5)	30% v/v 1,4-Dioxane	77.(G5)	
78.(G6)	0.001 M Magnesium chloride hexahydrate, 0.002 M Calcium chloride dihydrate	78.(G6)	0.05 M MOPS pH 7.0	78.(G6)	15% v/v (+/-)-2-Methyl-2,4-pentanediol	78.(G6)	
79.(G7)	0.08 M Strontium chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate	79.(G7)	0.04 M Sodium cacodylate trihydrate pH 7.0	79.(G7)	20% v/v (+/-)-2-Methyl-2,4-pentanediol	79.(G7)	0.012 M Spermine tetrahydrochloride
80.(G8)	0.08 M Sodium chloride	80.(G8)	0.04 M Sodium cacodylate trihydrate pH 7.0	80.(G8)	30% v/v (+/-)-2-Methyl-2,4-pentanediol	80.(G8)	0.012 M Spermine tetrahydrochloride
81.(G9)	0.04 M Lithium chloride, 0.08 M Strontium chloride hexahydrate	81.(G9)	0.04 M Sodium cacodylate trihydrate pH 7.0	81.(G9)	30% v/v (+/-)-2-Methyl-2,4-pentanediol	81.(G9)	0.012 M Spermine tetrahydrochloride
82.(G10)	0.04 M Lithium chloride, 0.08 M Strontium chloride hexahydrate, 0.02 M Magnesium chloride hexahydrate	82.(G10)	0.04 M Sodium cacodylate trihydrate pH 7.0	82.(G10)	30% v/v (+/-)-2-Methyl-2,4-pentanediol	82.(G10)	0.012 M Spermine tetrahydrochloride
83.(G11)	0.08 M Sodium chloride, 0.012 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	83.(G11)	0.04 M Sodium cacodylate trihydrate pH 7.0	83.(G11)	35% v/v (+/-)-2-Methyl-2,4-pentanediol	83.(G11)	0.012 M Spermine tetrahydrochloride
84.(G12)	0.012 M Sodium chloride, 0.08 M Potassium chloride	84.(G12)	0.04 M Sodium cacodylate trihydrate pH 7.0	84.(G12)	40% v/v (+/-)-2-Methyl-2,4-pentanediol	84.(G12)	0.012 M Spermine tetrahydrochloride
85.(H1)	0.08 M Sodium chloride, 0.02 M Barium chloride dihydrate	85.(H1)	0.04 M Sodium cacodylate trihydrate pH 7.0	85.(H1)	40% v/v (+/-)-2-Methyl-2,4-pentanediol	85.(H1)	0.012 M Spermine tetrahydrochloride
86.(H2)	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate	86.(H2)	0.04 M Sodium cacodylate trihydrate pH 7.0	86.(H2)	40% v/v (+/-)-2-Methyl-2,4-pentanediol	86.(H2)	0.012 M Spermine tetrahydrochloride
87.(H3)	0.08 M Potassium chloride, 0.02 M Barium chloride dihydrate	87.(H3)	0.04 M Sodium cacodylate trihydrate pH 7.0	87.(H3)	40% v/v (+/-)-2-Methyl-2,4-pentanediol	87.(H3)	0.012 M Spermine tetrahydrochloride
88.(H4)	0.08 M Potassium chloride, 0.02 M Magnesium chloride hexahydrate	88.(H4)	0.04 M Sodium cacodylate trihydrate pH 7.0	88.(H4)	50% v/v (+/-)-2-Methyl-2,4-pentanediol	88.(H4)	0.012 M Spermine tetrahydrochloride
89.(H5)	0.08 M Potassium chloride	89.(H5)	0.04 M Sodium cacodylate trihydrate pH 7.0	89.(H5)	60% v/v (+/-)-2-Methyl-2,4-pentanediol	89.(H5)	0.012 M Spermine tetrahydrochloride
90.(H6)	0.02 M Magnesium chloride hexahydrate, 0.002 M Cobalt(II) chloride hexahydrate	90.(H6)	0.05 M HEPES sodium pH 7.5	90.(H6)	2.0 M Ammonium sulfate	90.(H6)	0.001 M Spermine
91.(H7)	0.02 M Magnesium chloride hexahydrate	91.(H7)	0.05 M PIPES pH 7.5	91.(H7)	4% w/v Polyethylene glycol 8,000	91.(H7)	0.001 M Spermine
92.(H8)	0.015 M Magnesium chloride hexahydrate, 0.002 M Barium chloride dihydrate	92.(H8)	0.05 M PIPES pH 7.5	92.(H8)	7% v/v 2-Propanol	92.(H8)	0.0005 M Spermine
93.(H9)	0.02 M Magnesium chloride hexahydrate	93.(H9)	0.05 M PIPES pH 7.5	93.(H9)	10% w/v 1,6-Hexanediol	93.(H9)	0.001 M Spermine
94.(H10)	0.01 M Magnesium chloride hexahydrate	94.(H10)	0.05 M HEPES sodium pH 7.5	94.(H10)	15% v/v (+/-)-2-Methyl-2,4-pentanediol	94.(H10)	0.0015 M Spermine
95.(H11)	0.2 M Calcium chloride dihydrate	95.(H11)	0.05 M HEPES sodium pH 7.5	95.(H11)	28% v/v Polyethylene glycol 400	95.(H11)	0.002 M Spermine
96.(H12)	0.002 M Copper(II) chloride dihydrate	96.(H12)	0.05 M TRIS hydrochloride pH 8.5	96.(H12)	1.8 M Lithium sulfate monohydrate	96.(H12)	0.0005 M Spermine

◊ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

*Natrix HT contains ninety-six unique reagents. To determine the formulation of each reagent, simply read across the page.*



*Solutions for Crystal Growth*

34 Journey  
Aliso Viejo, CA 92656-3317 U.S.A.  
Tel: (949) 425-1321 • Fax: (949) 425-1011  
e-mail: tech@hmrll.com  
Website: www.hamptonresearch.com

© 1991 - 2018 Hampton Research Corp. all rights reserved  
Printed in the United States of America. This guide or parts thereof may not be reproduced in any form without the written permission of the publishers.