

# Updates to Detergent Screen HT & Detergent Screen Formulation



Detergent Screen HT and Detergent Screen

HR2-406 & HR2-408

## October 1, 2011

Reagent 57 (E9) in Detergent Screen (Detergent Screen HT) changed from 244.0 mM n-Octanoylsucrose to 250.0 mM SulfoBetaine 3-10.

Detergent Screen lots 240801, 240802, 240803 contain 244.0 mM n-Octanoylsucrose in reagent 57. Detergent lots 2408\*\* where \*\* is greater than or equal to 04 contain 250.0 mM SulfoBetaine 3-10.

Detergent Screen HT lots 240601, 240602, 240603, 240604, 240605 contain 244.0 mM n-Octanoylsucrose in reagent 57. Detergent lots 2406\*\* where \*\* is greater than or equal to 06 contain 250.0 mM SulfoBetaine 3-10.

## November 15, 2011

Reagent 2 (A2) in Detergent Screen (Detergent Screen HT) changed from 50.0 mM n-Dodecyl- $\beta$ -iminodipropionic acid, monosodium salt to 0.10 mM Lauryl maltose neopentyl glycol.

Detergent Screen lots 240801, 240802, 240803, 240804 contain 50.0 mM n-Dodecyl- $\beta$ -iminodipropionic acid, monosodium salt in reagent 2. Detergent Screen lots 2408\*\* where \*\* is greater than or equal to 05 contain 0.10 mM Lauryl maltose neopentyl glycol.

Detergent Screen HT lots 240601, 240602, 240603, 240604, 240605, 240606 contain 50.0 mM n-Dodecyl- $\beta$ -iminodipropionic acid, monosodium salt in reagent A2. Detergent Screen HT lots 2406\*\* where \*\* is greater than or equal to 07 contain 0.10 mM Lauryl maltose neopentyl glycol.

## November 2, 2012

Reagent 96 (H12) in Detergent Screen (Detergent Screen HT) changed from 70.0 mM LysoFos Choline 10 to 0.36 mM LysoFos Choline 14.

Detergent Screen lots 240801, 240802, 240803, 240804 contain 70.0 mM LysoFos Choline 10 in reagent 96. Detergent Screen lots 2408\*\* where \*\* is greater than or equal to 05 contain 0.36 mM LysoFos Choline 14.

Detergent Screen HT lots 240601, 240602, 240603, 240604, 240605, 240606, 240607 contain 70.0 mM LysoFos Choline 10 in reagent H12. Detergent Screen HT lots 2406\*\* where \*\* is greater than or equal to 08 contain 0.36 mM LysoFos Choline 14.

## January 14, 2013

Reagent Number 71 (F11) in Detergent Screen (Detergent Screen HT) changed from 6.6% w/v C-HEGA-8 changed to 45.0 mM Tripao.

Detergent Screen lots 240801, 240802, 240803, 240804, 240805 contain 70 mM LysoFos Choline 10 in reagent 71. Detergent Screen lots 2408\*\* where \*\* is greater than or equal to 06 contain 45.0 mM Tripao.

Detergent Screen HT lots 240601, 240602, 240603, 240604, 240605, 240606, 240607, 240608 contain 70 mM LysoFos Choline 10 in reagent F11. Detergent Screen HT lots 2406\*\* where \*\* is greater than or equal to 09 contain 45.0 mM Tripao.

Hampton Research  
34 Journey  
Aliso Viejo, CA 92656-3317 U.S.A.  
Tel: (949) 425-1321 • Fax: (949) 425-1611  
Technical Support e-mail: tech@hrmail.com  
Website: www.hamptonresearch.com

© 1991-2013 Hampton Research Corp. all rights reserved  
Printed in the United States of America. This guide or parts thereof may not be reproduced in any form without the written permission of the publishers.