




ALL DIMENSIONS IN MM



ITEM NO.	DESCRIPTION	INSTRUCTION	QTY.
No. 1	NON UVP Seal		1
No. 2	Protective Dust Cover (removeable)	Remove Dust Cover before imaging	1
No. 3	Pressure Adhesive	Activate adhesive by applying pressure	1
No. 4	Peelable Cover (remove first)	Remove seal from peelable cover	1



COPYRIGHT © 2020 SWISSCI LTD. ALL RIGHTS RESERVED



SWISSCI LTD  
 Lancaster Road, High Wycombe  
 Buckinghamshire, HP12 3PY  
 sales@swissci.com

TITLE HDP-PE-96T03 - Hanging Drop Seal NON UVP		PROJ NO.
DATE 10/08/2010	MATERIAL & CODE	LINE NO.
DRAWN BY AM	COLOUR & CODE	STATUS PRODUCTION
GENERAL TOLERANCE	SCALE 1 : 1	DRAWING NO. SHEET NO. 1

REV: 2020.A2.A