

Natrix™ 2 & Natrix HT™ Formulation Change

HAMPTON
RESEARCH

Solutions for Crystal Growth

HR2-117 & HR2-131

Formulation Change

Natrix 2 lots 211707 to 2117xx where xx > 07 and Natrix HT lots 213113-**-*** to 2131xx-**-*** where xx > 13 have the following formulation for Natrix 2 reagents 1, 2, 3, 4, and 26 (Natrix HT reagent E1, E2, E3, E4 and G2).

Hexamine cobalt(III) chloride concentration was lowered to 0.002 M from 0.020 M. Hampton Research received reports of possible Hexamine cobalt(III) chloride crystals in reagent and drop wells. No such reports had been received for reagents with 0.002 M Hexamine cobalt(III) chloride, so the concentration was lowered from 0.020 to 0.002 to reduce the likelihood for crystals of Hexamine cobalt(III) chloride.

HR2-917-01	Natrix 2 01 (E1)	0.04 M Lithium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 5.5, 30% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride
-------------------	-------------------------	--

HR2-917-02	Natrix 2 02 (E2)	0.08 M Sodium chloride, 0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 5.5, 35% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride
-------------------	-------------------------	---

HR2-917-03	Natrix 2 03 (E3)	0.012 M Sodium chloride, 0.08 M Potassium chloride, 0.04 M Sodium cacodylate trihydrate pH 5.5, 45% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride
-------------------	-------------------------	--

HR2-917-04	Natrix 2 04 (E4)	0.02 M Magnesium chloride hexahydrate, 0.04 M Sodium cacodylate trihydrate pH 5.5, 40% v/v (+/-)-2-Methyl-2,4-pentanediol, 0.002 M Hexamine cobalt(III) chloride
-------------------	-------------------------	--

HR2-917-26	Natrix 2 26 (G2)	0.02 M Magnesium chloride hexahydrate, 0.05 M MOPS pH 7.0, 55% v/v Tacsimate™ pH 7.0, 0.002 M Hexamine cobalt(III) chloride
-------------------	-------------------------	---

Hampton Research
34 Journey
Aliso Viejo, CA 92656-3317 U.S.A.
Tel: (949) 425-1321 • Fax: (949) 425-1611
Technical Support e-mail: tech@hrmail.com
Website: www.hamptonresearch.com

© 1991-2014 Hampton Research Corp. all rights reserved
Printed in the United States of America. This guide or
parts thereof may not be reproduced in any form without
the written permission of the publishers.